

## Product datasheet for **MR225357L4V**

### **Sv2b (NM\_153579) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Sv2b (NM_153579) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sv2b
Synonyms:	A830038F04Rik; AI850768; mKIAA0735
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153579
ORF Size:	2049 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225357).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_153579.4</a> , <a href="#">NP_705807.2</a>
RefSeq Size:	5548 bp
RefSeq ORF:	2052 bp
Locus ID:	64176
UniProt ID:	<a href="#">Q8BG39</a>
Cytogenetics:	7 D1
Gene Summary:	Probably plays a role in the control of regulated secretion in neural and endocrine cells. [UniProtKB/Swiss-Prot Function]



[View online »](#)