

Product datasheet for **MR225376**

G0s2 (NM_008059) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: G0s2 (NM_008059) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: G0s2
Synonyms: AI255151; AV006465
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR225376 ORF sequence
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAAGTGTGCAGGAGCTGATCCCTCTGGCCAAGGAGATGATGGCGCAGAAGCCCCGAGGGAAGCTAG
TGAAGCTATACGTGCTGGCAGTGTGCTGGCGCTCTCGGTGTGGTGTCTCGCCTAGTTGAGACGGTGTG
CAGCCCTTTCACAGCCGCCAGCCGTCTGCGCGACCAAGAGGCTGCAGTGGTGGAGCTGCGGGAAGCGTGT
GAACAGCAGTCCCTCCACAAGCAGGCCCTGCTGGCAGGAGCAAGGCACAGGAGGCGACCCTGTGCAGCC
GGGCCCTGTCCCTCCGGCAGCACGCCTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR225376 protein sequence
Red=Cloning site **Green**=Tags(s)
MESVQELIPLAKEMMAQKPRGKLVKLYVLGSQLALFGVVLGLVETVCSPTAASRLRDQEAAVVELREAC
EQQLHKQALLAGGKAQEATLCSRALSLRQHAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_008059

ORF Size: 312 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_008059.3](#), [NP_032085.1](#)

RefSeq Size: 847 bp

RefSeq ORF: 312 bp

Locus ID: 14373
UniProt ID: [Q61585](#)
Cytogenetics: 1 97.69 cM
MW: 11.1 kDa
Gene Summary: Promotes apoptosis by binding to BCL2, hence preventing the formation of protective BCL2-BAX heterodimers.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR225376