

Product datasheet for **MR225384L3V**

Limk1 (NM_010717) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Limk1 (NM_010717) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Limk1
Synonyms:	KIZ-1; Limk; LIMK-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010717
ORF Size:	1944 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225384).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010717.2 , NP_034847.1
RefSeq Size:	3323 bp
RefSeq ORF:	1944 bp
Locus ID:	16885
UniProt ID:	P53668
Cytogenetics:	5 G2



[View online »](#)

Gene Summary:

This gene encodes a member of the LIM kinase family of proteins. This protein is a serine/threonine kinase that regulates actin polymerization via phosphorylation and inactivation of the actin binding factor cofilin. This protein also stimulates axon growth and may play a role in brain development. Homozygous knockout mice for this gene exhibit reduced bone mass, abnormal neuronal morphology and altered synaptic function. Multiple pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015]