

## Product datasheet for **MR225386L3V**

### Ltbr (NM\_010736) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ltbr (NM_010736) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ltbr
Synonyms:	AI256028; CD18; Ltbr; LTbetaR; TNF-R-III; Tnfbr; TNFCR; TNFR-RP; TNFR2-RP; TNFRrp; Tnfrsf3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010736
ORF Size:	1245 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225386).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010736.3</a>
RefSeq Size:	2109 bp
RefSeq ORF:	1248 bp
Locus ID:	17000
UniProt ID:	<a href="#">P50284</a>
Cytogenetics:	6 59.32 cM
Gene Summary:	Receptor for the heterotrimeric lymphotoxin containing LTA and LTB, and for TNFS14/LIGHT. Promotes apoptosis via TRAF3 and TRAF5. May play a role in the development of lymphoid organs (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)