

Product datasheet for **MR225425**

Pak3 (NM_001195046) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pak3 (NM_001195046) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pak3
Synonyms:	PAK-3; Pak65alpha; Pak65beta; Stk4
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR225425 representing NM_001195046
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCTGACAGCTTGATAACGAAGAAAAACCCAGCTCCCCACTGAGGATGAACAGTAAACAACCGAG
 ACTCTTCAGCACTCAACACAGCTCCAAACCACTGCCCATGGCCCGGAAGAGAAGAATAAGAAAGCCAG
 GCTTCGCTCTATCTCCAGGAGGAGGGGATAAAACCAATAAGAAGAAAGAGAAAGCGCCAGAGATC
 TCTCTTCCTCAGACTTTGAGCATAAGTTCATGTGGGTTTTGATGCAGTCACCGGGGAATTCCTAACT
 CCCCTTCCAGACCTCTAGACCTGTGACGGTCGCTTCAAGTCAATCAGAGGGAAAAATGCCAGATCTCTA
 TGGCTCACAGATGTGCCAGGGAAGCTCCAGAGGGAATTCCTGAGCAATGGGCACGACTACTCCAAACC
 TCCAACATAACAAAGCTGGAACAGAAGAAGAACCCACAAGCTGTTCTTGATGTTCTCAAGTTCTACGACT
 CCAAAGAAACGGTCAACAACAGAAATATATGAGCTTTACGTGAGGAGATAAAAGTGCCCATGGATACAT
 AGCAGCACATCAGTCGAATACCAAAACAGCTTCGGAACCTCCTTTGGCTCCTCTGTATCTGAAGAAGAG
 GATGAAGAAGAGGAAGAGGAAGAAGATGACAATGAACCTCCGCCTGTATTGCACCAAGACCAGAGCATA
 CAAAATCAATCTATACTCGCTCCGTGGTTGAGTCAATTGCTTCACCAGCTGCACCAATAAAGAAGATAT
 CCCACCTTCTGCTGAGAAATGCCAATCCACCACCTTTGTACAGGAATACAGATCGGCAAGAAAAAAGTCC
 AAGATGACGGATGAGGAGATCCTAGAGAAGCTGAGAAGCATTGTGAGTGTGGGGACCCAAAGAAGAAAT
 ATACGAGATTTGAAAAAATGGCCAAGGGGCATCAGGAAGTGTACACAGCACTAGACATTGCGACTGG
 ACAAGAGGTGGCCATAAAGCAATGAACCTTCAACAGCAGCCAAAAAGGAATTAATTATTAATGAAATT
 CTTGTCTAGGGGAAAAAAGAACCCTTCAATTGCAATTATTTAGATAGCTACTAGTGGGTGATGAAC
 TGTGGGTAGTCATGGAATACTTGGCTGGTGGCTCTTTGACTGATGTGGTCACAGAAACCTGTATGGATGA
 AGGACAGATAGCAGCTGTCTGTAGAGAGTGCCTCCAAGCTTTGGATTTCTTGCACTCAAACCAAGTGATC
 CACAGAGACATAAAGAGCGACAACATTCTCCTCGGGATGGATGGTTCTGTTAAACTGACTGATTTTGGGT
 TCTGTGCTCAAATCACTCCTGAGCAAAGTAAACGAAGCACTATGGTGGGAACTCCCTATTGGATGGCACC
 TGAAGTGGTAACTCGAAAAGCTTATGGTCCAAAAGTTGATATCTGGTCTCTGGGAATCATGGCCATTGAA
 ATGGTGGAAAGTGAACCCCTTACCTTAATGAAAATCCACTCAGGGCCTTATATCTGATAGCCACTAATG
 GGACCCAGAGCTCCAGAATCCTGAGAGACTGTCAGCTGTATTCCGTGACTTCTAAATCGCTGTCTTGA
 GATGGATGTGGATAGAAGAGGATCTGCTAAGGAGCTTTTGCAGCATCCGTTTTTAAATAGCCAAGCCT
 CTGTCTAGTCTACTCCTCTGATTATCGCCGAAAGGAAGCAATTAAGAACAGTAGCCGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR225425 representing NM_001195046
 Red=Cloning site Green=Tags(s)


MSDSLDNEEKPPAPPLRMNSNDRSSALNHSSKPLPMAPEEKKKARLRSIFPGGGDKTNKKKEKERPEI
 SLPSPDFEHTIHVGFDAVTGEFTNSPFQTSRPVTVASSQSEGKMPDLYGSQMCPGKLEPIPEQWARLLQT
 SNITKLEQKKNPQAVLDVLFYDSKETVNNQKYSFTSGDKSAHGYYIAAHQSNTKTASEPPLAPPVSEEE
 DEEEEEEDDNEPPPVIAPRPEHTKSIYTRSVVESIASPAAPNKEDIPPSAENANSTTL YRNTDRQRKKS
 KMTDEEILEKLRISIVSGDPKKKYTRFEKIGQGASGTVYALDIATGQEVAIKQMNLLQQPKKELINEI
 LVMRENKPNINYLDSYLVGDELWVMEYLAGGSLTDVVTETCMDEGQIAAVCRECLQALDFLHSNQVI
 HRDIKSDNILLGMDGSVKLTDGFCAQITPEQSKRSTMVGTPTYWMAPEVVTRKAYGPKVDIWSLGIMAIE
 MVEGEPYLNENPLRALYLIATNGTPELQNPRLSAVFRDFLNRCLMDVDRRGSAKELLQHPFLKLAKP
 LSSLTPLIIAAKEAIKNSSR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001195046

ORF Size: 1740 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001195046.1 , NP_001181975.1
RefSeq Size:	8242 bp
RefSeq ORF:	1743 bp
Locus ID:	18481
Cytogenetics:	X F2
MW:	65.1 kDa
Gene Summary:	<p>Serine/threonine protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell migration, or cell cycle regulation. Plays a role in dendrite spine morphogenesis as well as synapse formation and plasticity (PubMed:25851601). Acts as downstream effector of the small GTPases CDC42 and RAC1. Activation by the binding of active CDC42 and RAC1 results in a conformational change and a subsequent autophosphorylation on several serine and/or threonine residues. Phosphorylates MAPK4 and MAPK6 and activates the downstream target MAPKAPK5, a regulator of F-actin polymerization and cell migration. Additionally, phosphorylates TNNI3/troponin I to modulate calcium sensitivity and relaxation kinetics of thin myofilaments. May also be involved in early neuronal development.[UniProtKB/Swiss-Prot Function]</p>