

Product datasheet for **MR225439L4V**

Zswim2 (NM_027964) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zswim2 (NM_027964) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zswim2
Synonyms:	1700025P14Rik; 4933437F18Rik; MEX
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027964
ORF Size:	1893 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225439).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027964.2 , NP_082240.1
RefSeq Size:	2062 bp
RefSeq ORF:	1896 bp
Locus ID:	71861
UniProt ID:	Q9D9X6
Cytogenetics:	2 D
Gene Summary:	E3 ubiquitin-protein ligase involved in the regulation of Fas-, DR3- and DR4-mediated apoptosis. Functions in conjunction with the UBE2D1, UBE2D3 and UBE2E1 E2 ubiquitin-conjugating enzymes.[UniProtKB/Swiss-Prot Function]


[View online »](#)