

Product datasheet for **MR225476L3V**

Crabp1 (NM_013496) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Crabp1 (NM_013496) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Crabp1
Synonyms:	AI326249; Crabp-1; CRABP-I; Crabpl; Rbp-5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013496
ORF Size:	411 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225476).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013496.2 , NP_038524.1
RefSeq Size:	741 bp
RefSeq ORF:	414 bp
Locus ID:	12903
UniProt ID:	P62965
Cytogenetics:	9 29.76 cM
Gene Summary:	Cytosolic CRABPs may regulate the access of retinoic acid to the nuclear retinoic acid receptors.[UniProtKB/Swiss-Prot Function]



[View online »](#)