

## Product datasheet for **MR225485L1V**

### **Pmp22 (NM\_008885) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pmp22 (NM_008885) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pmp22
Synonyms:	22kDa; Gas-3; HNPP; Tr; trembler
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_008885
ORF Size:	480 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225485).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008885.2</a> , <a href="#">NP_032911.1</a>
RefSeq Size:	1801 bp
RefSeq ORF:	483 bp
Locus ID:	18858
UniProt ID:	<a href="#">P16646</a>
Cytogenetics:	11 38.99 cM
Gene Summary:	Might be involved in growth regulation, and in myelination in the peripheral nervous system.[UniProtKB/Swiss-Prot Function]



[View online »](#)