

Product datasheet for **MR225501**

Fzd3 (NM_021458) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fzd3 (NM_021458) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fzd3
Synonyms:	AU020229; D930050A07Rik; Fz3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR225501 representing NM_021458
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGCTGTGAGCTGGATTGTCTTTGATCTTTGGCTCTTACTGTGTTTCTGGGGCAGATAGTGGGCACA
GTTTGTCTTTCTTGGAACCTATAACCTTGAGGATGTGCCAAGATTGCCTTACAATACTACCTTCATGCC
TAATCTTCTGAACCATTATGACCAACAGACTGCAGCTTTAGCAATGGAGCCCTTCCACCCTATGGTGAAC
CTGGATTGTTCTCGGATTTTCGGCCATTTCTTTGTGCACTCTATGCCCTATTTGTATGGAATATGGAC
GTGTCACACTTCCCTGCCGTAGGCTGTGTGAGCGTGCCTATAGCGAGTGTCAAACCTCATGGAGATGTT
TGGTGTCCCGTGGCCTGAAGATATGGAGTGCAGTAGGTTCCAGATTGTGATGAGCCATATCCCGACTT
GTGGATTTGAATTTAGTTGGAGATCCAACCTGAAGGAGCCCAAGTGCAGTGCAGAGGGACTATGGTTTTT
GGTGTCCAGAGAGTTAAAAATGATCCTGATCTTGGCTATTCTTTCTGCACGTGCGAGATTGTCGCC
ACCATGTCCCAATATGTACTTCAGGAGAGAAGAACTGTCAATTTGCTCGCTATTTTCATAGGCTGATTTCA
ATCATTGCTCTCTGCCACATTGTTTACTTTTTAACCTTTCTAATTGACGTGACAAGATTCCGTTACC
CTGAAAGACCTATCATATTTTATGCAGTCTGCTACATGATGGTGTCAATATTTTCTTATTGGGTTTTT
GCTGGAGGACCGAGTAGCCTGCAATGCATCTAGCCCTGCACAGTATAAGGCTTCTACAGTGACACAAGGA
TCTCACAATAAGGCTGTACCATGCTCTTTATGGTACTATATTTTTTCACTATGGCTGGCAGTGTATGGT
GGGTAATCTTACCATCATATGGTTTTTAGCAGCTGTGCCAAAGTGGGGCAGTGAAGCTATTGAGAAGAA
AGCATTGCTGTTTATGCCAGTGCCTGGGGCATCCCCGGAACCTAACTATCATCTTTTAGCGATGAAT
AAAATTGAAGGTGACAATATTAGTGGCGTGTGTTTTGTCGGCCTACGACGTTGATGCATTAAGATATT
TCGTTCTCGCTCCCCTCTGCCTGTATGTGGTAGTTGGGTTTTCTCCTTTTAGCCGGCATTATATCCCT
AAACAGAGTTCGGATTGAGATCCCATTAGAAAAGGAAAACCAAGATAAGTTAGTGAAGTTCATGATCCGG
ATTGGTGTTCAGCATTCTCTACCTTGTGCCACTCTTGGTTGTAATTGGATGTTACTTTTATGAGCAAG
CTTACCGCGCATCTGGGAGACAACATGGATCCAGGAACGCTGCAGAGAGTATCACATTCATGTCGGTA
CCAGGTTACTCAGATGAGTCGTCCAGACCTGATTCTTTCTGATGAAGTATCTCATGGCTCTCATAGTT
GGGATTCCTCTATATTTGGGTTGGAAGCAAAAAGACATGCTTTGAATGGGCCAGTTTTTCCATGGGC
GTAGGAAAAAAGAGATAGTGAATGAGAGCCGGCAGGTGCTCCAGGAACCTGACTTTGCTCAGTCACTCCT
GAGGGACCCAAATACTCCAATTATAAGAAAATCAAGAGGAACCTCAGTCAAGGGACATCCACACATGCT
TCTTCAACTCAGCTGGCCATGGTGGATGACCAAGAAGCAAAGCAGGGAGTGTCCACAGCAAAGTGAGCA
GCTACCATGGCAGCCTCCACAGGTACGGGATGGCAGGTACACTCCCTGCAGTTACCGAGGAATGGAGGA
GAGACTACCTCACGGCAGCATGTCACGGCTGACGGATCATTCCAGGCACAGTAGTTCTCATCGGCTCAAC
GAGCAGTCCCACACAGCAGCATCCGAGACCTCAGTAACAACCCCATGACTCACATTACACATGGCACCA
GCATGAACCGTGTATTGAGGAGGATGGAACCAAGTGT

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR225501 representing NM_021458
 Red=Cloning site Green=Tags(s)


MAVSWIVFDLWLLTVFLGQIGGHSLSFCEPITLRMCQDLPYNTTFMPNLLNHYDQQTAAALAMEPFHPMVN
 LDCSRDFRPFLCALYAPICMEYGRVTLPCRRLCQRAYSECSKLMEMFGVPWPEDMECSRFPDCDEPYPRL
 VDLNLVGDPTGAPVAVQRDYGFWCPRELKIDPDLGYSFLHVRDCSPPCPNMYFRREELSFARYFIGLIS
 IICLSATLFTFLTFLIDVTRFRYPERPIIFYAVCYMMVSLIFFIGFLEDRVACNASSPAQYKASTVTQG
 SHNKACTMLFMVLYFFTMAGSVVWVILTIWFLAAVPKWGSEAIEKKALLFHASAWGIPGTLTIILLAMN
 KIEGDNISGVCFVGLYDVDALRYFVLAFLCLYVVVGVSLLLAGIISLNRVRIEIPLEKENQDKLVKFMIR
 IGVFSILYLVPLLVVIGCYFYEQAIRGIWETTWIQERCREYHPCPYQVTQMSRPDLILFLMKYLMALIV
 GIPSIFWVGSKTCFEWASFFHGRRKKEIVNESRQVLQEPDFAQSLLRDPNTPIIIRKSRGTSTQGTSTHA
 SSTQLAMVDDQSKAGSVHSKSVSSYHGSLHRSRDGRYTPCSYRGMEEERLPHGMSRLTDHSRHSSSHRLN
 EQSRHSSIRDLSNPMTHITHGTSMNRVIEEDG TSA

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1164_e10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_021458

ORF Size: 1998 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021458.2](#), [NP_067433.1](#)

RefSeq Size: 12742 bp

RefSeq ORF: 2001 bp

Locus ID: 14365


UniProt ID: [Q61086](#)

Cytogenetics: 14 34.09 cM

MW: 76.7 kDa

Gene Summary:

Receptor for Wnt proteins. Most of frizzled receptors are coupled to the beta-catenin canonical signaling pathway, which leads to the activation of disheveled proteins, inhibition of GSK-3 kinase, nuclear accumulation of beta-catenin and activation of Wnt target genes. A second signaling pathway involving PKC and calcium fluxes has been seen for some family members, but it is not yet clear if it represents a distinct pathway or if it can be integrated in the canonical pathway, as PKC seems to be required for Wnt-mediated inactivation of GSK-3 kinase. Both pathways seem to involve interactions with G-proteins. Activation by Wnt5A stimulates PKC activity via a G-protein-dependent mechanism. Involved in transduction and intercellular transmission of polarity information during tissue morphogenesis and/or in differentiated tissues. Plays a role in controlling early axon growth and guidance processes necessary for the formation of a subset of central and peripheral major fiber tracts. Required for the development of major fiber tracts in the central nervous system, including: the anterior commissure, the corpus callosum, the thalamocortical, corticothalamic and nigrostriatal tracts, the corticospinal tract, the fasciculus retroflexus, the mammillothalamic tract, the medial lemniscus, and ascending fiber tracts from the spinal cord to the brain. In the peripheral nervous system, controls axon growth in distinct populations of cranial and spinal motor neurons, including the facial branchiomotor nerve, the hypoglossal nerve, the phrenic nerve, and motor nerves innervating dorsal limbs. Involved in the migration of cranial neural crest cells. May also be implicated in the transmission of sensory information from the trunk and limbs to the brain. Controls commissural sensory axons guidance after midline crossing along the anterior-posterior axis in the developing spinal cord in a Wnt-dependent signaling pathway. Together with FZD6, is involved in the neural tube closure and plays a role in the regulation of the establishment of planar cell polarity (PCP), particularly in the orientation of asymmetric bundles of stereocilia on the apical faces of a subset of auditory and vestibular sensory cells located in the inner ear. Promotes neurogenesis by maintaining sympathetic neuroblasts within the cell cycle in a beta-catenin-dependent manner. [UniProtKB/Swiss-Prot Function]

Product images:

Circular map for MR225501