

## Product datasheet for **MR225518L4V**

### Kctd8 (NM\_175519) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Kctd8 (NM_175519) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kctd8
Synonyms:	A730087N02Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175519
ORF Size:	1431 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225518).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_175519.5</a> , <a href="#">NP_780728.3</a>
RefSeq Size:	2889 bp
RefSeq ORF:	1431 bp
Locus ID:	243043
UniProt ID:	<a href="#">Q50H33</a>
Cytogenetics:	5 C3.1


[View online »](#)

**Gene Summary:**

Auxiliary subunit of GABA-B receptors that determine the pharmacology and kinetics of the receptor response. Increases agonist potency and markedly alter the G-protein signaling of the receptors by accelerating onset and promoting desensitization.[UniProtKB/Swiss-Prot Function]