

Product datasheet for **MR225519L3V**

Shc3 (NM_009167) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Shc3 (NM_009167) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Shc3
Synonyms:	N-Shc; Rai; ShcC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009167
ORF Size:	1422 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225519).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009167.3
RefSeq Size:	1540 bp
RefSeq ORF:	1425 bp
Locus ID:	20418
UniProt ID:	Q61120
Cytogenetics:	13 26.12 cM
Gene Summary:	Signaling adapter that couples activated growth factor receptors to signaling pathway in neurons. Involved in the signal transduction pathways of neurotrophin-activated Trk receptors in cortical neurons (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)