

Product datasheet for **MR22546L3V**

Anxa7 (NM_001110794) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Anxa7 (NM_001110794) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Anxa7
Synonyms:	AI265384; AI316497; Anx7; synexin
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001110794
ORF Size:	1392 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR22546).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001110794.1 , NP_001104264.1
RefSeq Size:	2850 bp
RefSeq ORF:	1392 bp
Locus ID:	11750
UniProt ID:	Q07076
Cytogenetics:	14 11.53 cM
Gene Summary:	Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.[UniProtKB/Swiss-Prot Function]



[View online »](#)