

## Product datasheet for **MR225549L4V**

### **Ptger1 (NM\_013641) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ptger1 (NM_013641) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ptger1
Synonyms:	42kDa; EP1; Ptgerep1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013641
ORF Size:	1215 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225549).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013641.2</a>
RefSeq Size:	2082 bp
RefSeq ORF:	1218 bp
Locus ID:	19216
UniProt ID:	<a href="#">P35375</a>
Cytogenetics:	8 40.22 cM



[View online »](#)

**Gene Summary:**

Receptor for prostaglandin E2 (PGE2). The activity of this receptor is mediated by G(q) proteins which activate a phosphatidylinositol-calcium second messenger system. May play a role as an important modulator of renal function. Implicated the smooth muscle contractile response to PGE2 in various tissues.[UniProtKB/Swiss-Prot Function]