

Product datasheet for **MR225550L3V**

Gabrg2 (NM_177408) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gabrg2 (NM_177408) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gabrg2
Synonyms:	GABAA-R; Gabrg-2; gamma2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177408
ORF Size:	1398 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225550).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177408.6
RefSeq Size:	3511 bp
RefSeq ORF:	1401 bp
Locus ID:	14406
UniProt ID:	P22723
Cytogenetics:	11 24.8 cM



[View online »](#)

Gene Summary:

This gene encodes a gamma-aminobutyric acid (GABA)-A receptor subunit, which is a member of the ligand-gated ion channel family. GABA is the major inhibitory neurotransmitter in the adult central nervous system, and conversely exhibits an excitatory function during development. GABA-A receptors are pentameric, consisting of proteins from several subunit classes: alpha, beta, gamma, delta and rho. This gene encodes one of three gamma subunits in mammals, which contain the binding site for benzodiazepine drugs. Several mutations in this gene are associated with epileptic seizures, and genetic knockdown is associated with anxiety behavior. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013]