

Product datasheet for MR225565

Irf1 (NM_001159393) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Irf1 (NM_001159393) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Irf1
Synonyms:	AU020929; Irf-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR225565 representing NM_001159393 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCAATCACTCGAATGCGGATGAGACCTGGCTAGAGATGCAGATTAATTCACCAAAATCCCAGGGC
TGATCTGGATCAATAAAGAAGAGATGATCTTCCAGATTCCATGGAAGCACGCTGCTAAGCACGGCTGGGA
CATCAACAAGGATGCCTGTCTGTTCCGGAGCTGGGCCATTCACACAGGCCGATACAAAGCAGGAGAAAAA
GAGCCAGATCCCAAGACATGGAAGGCAAACCTCCGTTGTGCCATGAACCTCCGAGACATCGAGGAAG
TGAAGGATCAGAGTAGGAACAAGGCAGCTCTGCTGTGCGGGTGTACCGGATGCTGCCACCCCTCACCAG
GAACCAGAGGAAAAGAGAGAAAGTCCAAGTCCAGCCGAGACACTAAGAGCAAACCAAGAGGAAGCTGTGT
GGAGATGTTAGCCCGGACACTTTCTCTGATGGACTCAGCAGCTCTACCCTACCTGATGACCACAGCAGTT
ACACCACTCAGGGCTACCTGGGTCAGGACTTGGATATGAAAAGGGACATAACTCCAGCACTGTCACCGTG
TGTCGTCAGCAGCAGTCTCTCTGAGTGGCATATGCAGATGGACATTATACCAGATAGCACCAGTCTG
TATAACCTACAGGTGTCACCCATGCCTTCCACCTCCGAAGCCGCAACAGACGAGGATGAGGAAGGGAAGA
TAGCCGAAGACCTTATGAAGCTCTTGAACAGTCTGAGTGGCAGCCGACACATCGATGGCAAGGGATA
CTTGCTCAATGAGCCAGGGACCCAGCTCTCTTCTGTCTATGGAGACTTCAGCTGCAAAGAGGAACCCAGAG
ATTGACAGCCCTCGAGGAAACCTCCTGATGGGAGTCTTTTGCTGGCTTTCTGCTGGGCTTCAGCTGAGC
AT


ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001159393.1, NP_001152865.1</u>
RefSeq Size:	2128 bp
RefSeq ORF:	915 bp
Locus ID:	16362
Cytogenetics:	11 32.0 cM
MW:	35.1 kDa
Gene Summary:	<p>Transcriptional regulator which displays a remarkable functional diversity in the regulation of cellular responses. These include the regulation of IFN and IFN-inducible genes, host response to viral and bacterial infections, regulation of many genes expressed during hematopoiesis, inflammation, immune responses and cell proliferation and differentiation, regulation of the cell cycle and induction of growth arrest and programmed cell death following DNA damage. Stimulates both innate and acquired immune responses through the activation of specific target genes and can act as a transcriptional activator and repressor regulating target genes by binding to an interferon-stimulated response element (ISRE) in their promoters. Its target genes for transcriptional activation activity are: genes involved in anti-viral response, such as IFN-alpha/beta, DDX58/RIG-I, TNFSF10/TRAIL, OAS1/2, PIAS1/GBP, EIF2AK2/PKR and RSAD2/viperin; antibacterial response, such as NOS2/INOS; anti-proliferative response, such as p53/TP53, LOX and CDKN1A; apoptosis, such as BBC3/PUMA, CASP1, CASP7 and CASP8; immune response, such as IL7, IL12A/B and IL15, PTGS2/COX2 and CYBB; DNA damage responses and DNA repair, such as POLQ/POLH; MHC class I expression, such as TAP1, PSMB9/LMP2, PSME1/PA28A, PSME2/PA28B and B2M and MHC class II expression, such as CIITA. Represses genes involved in anti-proliferative response, such as BIRC5/survivin, CCNB1, CCNE1, CDK1, CDK2 and CDK4 and in immune response, such as FOXP3, IL4, ANXA2 and TLR4. Stimulates p53/TP53-dependent transcription through enhanced recruitment of EP300 leading to increased acetylation of p53/TP53. Plays an important role in immune response directly affecting NK maturation and activity, macrophage production of IL12, Th1 development and maturation of CD8+ T-cells. Also implicated in the differentiation and maturation of dendritic cells and in the suppression of regulatory T (Treg) cells development. Acts as a tumor suppressor and plays a role not only in antagonism of tumor cell growth but also in stimulating an immune response against tumor cells.[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR225565