

Product datasheet for **MR225583L4V**

Gria1 (NM_001113325) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gria1 (NM_001113325) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gria1
Synonyms:	2900051M01Rik; Glr-1; Glr1; GluA1; Glur-1; GluR-A; gluR-K1; Glur1; GluRA; HIPA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001113325
ORF Size:	2721 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225583).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001113325.2 , NP_001106796.1
RefSeq Size:	5755 bp
RefSeq ORF:	2724 bp
Locus ID:	14799
UniProt ID:	P23818
Cytogenetics:	11 34.51 cM



[View online »](#)

Gene Summary:

Ionotropic glutamate receptor. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous system. Binding of the excitatory neurotransmitter L-glutamate induces a conformation change, leading to the opening of the cation channel, and thereby converts the chemical signal to an electrical impulse. The receptor then desensitizes rapidly and enters a transient inactive state, characterized by the presence of bound agonist. In the presence of CACNG4 or CACNG7 or CACNG8, shows resensitization which is characterized by a delayed accumulation of current flux upon continued application of glutamate.[UniProtKB/Swiss-Prot Function]