

## Product datasheet for **MR225584L3V**

### Sell (NM\_001164059) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sell (NM_001164059) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sell
Synonyms:	AI528707; CD62L; L-selectin; LECAM-1; Lnhr; Ly-22; Ly-m22; Lyam-1; Lyam1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164059
ORF Size:	1008 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225584).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001164059.1</a> , <a href="#">NP_001157531.1</a>
RefSeq Size:	2213 bp
RefSeq ORF:	1011 bp
Locus ID:	20343
Cytogenetics:	1 71.37 cM
Gene Summary:	Calcium-dependent lectin that mediates cell adhesion by binding to glycoproteins on neighboring cells. Mediates the adherence of lymphocytes to endothelial cells of high endothelial venules in peripheral lymph nodes (PubMed:1693096). Promotes initial tethering and rolling of leukocytes in endothelia (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)