

Product datasheet for **MR225593L1V**

Nid1 (NM_010917) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nid1 (NM_010917) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nid1
Synonyms:	A630025O17; entactin; entactin-1; Nid; nidogen-1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_010917
ORF Size:	3738 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225593).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010917.2 , NP_035047.2
RefSeq Size:	6046 bp
RefSeq ORF:	3738 bp
Locus ID:	18073
UniProt ID:	P10493
Cytogenetics:	13 5.26 cM
Gene Summary:	Sulfated glycoprotein widely distributed in basement membranes and tightly associated with laminin. Also binds to collagen IV and perlecan. It probably has a role in cell-extracellular matrix interactions.[UniProtKB/Swiss-Prot Function]


[View online »](#)