

Product datasheet for **MR225611A1V**

Mouse Oxt (NM_011025) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Oxt (NM_011025) AAV Particle
Tag: Myc-DDK
Symbol: Oxt
Synonyms: OT; Ox; Oxy
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR225611 representing NM_011025
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCCTGCCCCAGTCTCGCTTGCTGCCTGCTTGCTTACTGGCTCTGACCTCGGCCTGCTACATCCAGA
 ACTGCCCCCTGGGCGGCAAGAGGGCTGTGCTGGACCTGGATATGCGCAAGTGTCTCCCTGCGGCCCGGG
 CGGCAAGGACGCTGCTTCGGACCAAGCATCTGCTGCGCGGACGAGCTGGGTGCTTCGTGGGCACCGCC
 GAGGCGCTGCGCTGCCAGGAGGAGAATACCTGCCTTCGCCCTGCCAGTCTGGCCAGAAGCCCTGCGGGA
 CGGAGGCCGCTGCGCCGCCACAGGCATCTGCTGCAGCCCGGATGGCTGCCGCACGGACCCCGCCTGCGA
 CCCTGAGTCTGCCTTCTCGGAGCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR225611 representing NM_011025
 Red=Cloning site Green=Tags(s)

MACPSLACLLGLLALTSACYIQNCPLGGKRAVLDLDMRKCLPCPGGKGRCFGPSICCADELGCFVGT
 EALRCQEENYLPSPCQSGQKPCGSGGRCAATGICCPDGCRTDPACDPESAFSER

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_011025



[View online »](#)

ORF Size:	375 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_011025.3</u> , <u>NP_035155.1</u>
RefSeq Size:	471 bp
RefSeq ORF:	378 bp
Locus ID:	18429
UniProt ID:	<u>P35454</u>
Cytogenetics:	2 63.24 cM
MW:	13.3 kDa