

Product datasheet for **MR225685**

Chat (NM_009891) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Chat (NM_009891) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Chat
Synonyms:	B230380D24Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR225685 representing NM_009891
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCTATCCTGGAAAAGTCCCCCAAAGATGCCTGTACAAGCTTCTAGCTGTGAGGAGGTGCTGGACT
TACCTAAGTTGCCAGTGCCCCACTGCAGCAAACCTGGCCACCTACCTTCAGTGCATGCAACACCTGGT
ACCTGAAGAGCAGTTCAGGAAGAGCCAGGCCATTGTGAAGCGTTTGGGGCCCTGGTGGCCTGGGTGAG
ACCCTGCAGGAAAAGCTCTTGGAGAGACAGGAGAAGACAGCCAATTGGGTCTCTGAATACTGGTGAATG
ACATGTATCTAAACAACCGCCTGGCCCTGCCAGTCAACTCTAGCCCTGCTGTGATCTTTGCTCGGCAGCA
CTTCCAAGACACCAATGACCAGCTAAGGTTTGCAGCCAGCCTCATCTCTGGTGTGCTTAGCTACAAGGCT
CTGCTGGACAGCAATCCATTCCACTGACTGGCCAAGGGCAGCTCTCAGGGCAGCCTCTCTGTATGA
AGCAGTACTACAGACTCTTCTCATCATACCGGCTTCTGGCCATACCCAGGACACACTGGTGGCCAGAA
GAGCAGTATCATGCCTGAGCCCGAGCATGTCATCGTGGCCTGCTGCAACCAGTTCCTTTGCTTGGATGTT
GTCAATTAATTCGCCGCTCAGTGAGGGTGATCTGTTCACTCAGTTGAGAAAGATAGTCAAATGGCGT
CCAACGAGGATGAACGCCTGCCTCCAATCGGCCTGCTGACGTCAGACGGGAGGAGCGAGTGGGCCAAGGC
CAGGACGGTCTCTTAAAGACTCCACCAACCGGGACTCCCTGGACATGATCGAGCGCTGCATCTGCCTG
GTATGCCTGGATGGTCCAGGCACTGGAGACCTCAGTGACACACAGGGCCCTCCAGTCCCTCATGGTG
GAGGCTGCAGCTTGAATGGAGCGAATCGTTGGTATGACAAGTCCCTGCAGTTTGTGGTGGGCCGAGATGG
CACCTGCGGTGTGGTGTGTGAGCACTCCCTTTTATGGCATCGTCTGGTGCAGTGCACGGAGCACCTG
CTGAAACATATGATGACAGGCAACAAGAAGCTCGTCCGAGCTGACTCAGTGAGTGAATCCCTGCTCCCA
GAAGGCTGAGGTGAAATGTTCCCAAGAACTCAAGCCATCTCGCCTCCTCGGCAGAGAACTTCAAAG
AATCGTAAAGAATCTGGATTTCAATGTTTATAAGTTTGACAACATGAGGAAAACATTTATTAAGAAGCAG
AAATGCAGTCTGATGGCTTCATCCAGGTGGCCCTCCAGCTGGCTTACTACAGGCTTTACCAGAGGCTGG
TGCCACCTATGAGAGTGCATCCATCCGCCGCTTCCAGGAAGGTCGGGTGGACAACATCAGATCGGCCAC
TCCTGAGGCTCTGGCTTTTGTGCAAGCCATGACTGACCACAAGGCTGCCGTGCTGGCTTCTGAGAACTG
CAGCTGCTGCAGAGGGCCATCCAGGCCAAACTGAGTACACAGTCATGGCCATAACCGGCATGGCCATTG
ACAACCATCTTCTGGCACTGAGGGAGCTGGCCCGAGACCTGTGCAAAGAGCCACCTGAGATGTTTATGGA
TGAACATACCTGATGAGCAACCGGTTTATTCTCTCCACCAGCCAGGTGCCACGACCATGGAGATGTTT
TGCTGTTATGGCCCTGTGGTACCAATGGGTATGGAGCCTGTTATAACCCCGAGCCTGAGGCCATCACCT
TCTGCATCTCCAGCTTTCACGGCTGCAAAGAGACCTCATCTGTGGAGTTTGCAGAAGCGGTGGGAGCGAG
CCTTGTGACATGAGAGCCTCTGTAGTTCGAGGCAGCCTGCTGAAGGCAAGCCACCAACAGCAAAGGAA
AGAGCTAGAGGCCCAACCAAGCAAGCAATCTTGACTACTCCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR225685 representing NM_009891
 Red=Cloning site Green=Tags(s)

```


MPILEKVPKMPVQASSCEEVLDLPKLPVPLQQTALATYLQCMQHLVPPEQFRKSQAIVKRFGAPGGLGE
TLQEKLLEKQKQKTANVWSEYWLNDMYLNNRLALPVNSSPAVIFARQHFDQTNQDLRFAASLISGVL SYKA
LLDSQSIPTDWAKGQLSGQPLCMKQYYRLFSSYRLPGHTQDTLVAQKSSIMPEPEHVIVACCNQFFVLDV
VINFRRLSEGDLFTQLRKIVKMASNEDERLPPIGLLTSDGRSEWAKARTVLLKDDSTNRDSLMIERCICL
VCLDGPPTGDLSDTHRALQLLHGGGCSLNGANRWYDKSLQFVVGRDGTGCVVCEHSPFDGIVLVQCTEHL
LKHMMTGNKKLVRAHSVSELPAPRRLRWKCSPETQGHLLASSAEKLQRIKVNLDLFIYKFDNYGKTFIKKQ
KCSPDGF IQVALQLAYYRLYQRLVPTYESASIRRFQEGRVDNIRSATPEALAFVQAMTDHKA AVLASEKL
QLLQRAIQAQTEYTVMAITGMAIDNHLLALRELARDLCKEPPMEMFDETYLMSNRFILSTSQVPTTMMEMF
CCYGPVVPNGYGACYNPQPEAITFCISSFHGCKETSSVEFAEAVGASLVDMRDLCSSRQPAEGKPPTAKE
RARGPTKPSNLDYSH
  
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9075_a04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_009891

ORF Size: 1935 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009891.2](#), [NP_034021.1](#)

RefSeq Size: 2810 bp


RefSeq ORF: 1938 bp

Locus ID: 12647

Cytogenetics: 14 19.4 cM

MW: 72.3 kDa

Product images:



Circular map for MR225685