

Product datasheet for **MR225686L4V**

Lrat (NM_023624) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrat (NM_023624) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrat
Synonyms:	1300010A18Rik; AI449251
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_023624
ORF Size:	696 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225686).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023624.4 , NP_076113.1
RefSeq Size:	5349 bp
RefSeq ORF:	696 bp
Locus ID:	79235
UniProt ID:	Q9J160
Cytogenetics:	3 E3



[View online »](#)

Gene Summary:

Transfers the acyl group from the sn-1 position of phosphatidylcholine to all-trans retinol, producing all-trans retinyl esters (By similarity). Retinyl esters are storage forms of vitamin A (By similarity). LRAT plays a critical role in vision (By similarity). It provides the all-trans retinyl ester substrates for the isomerohydrolase which processes the esters into 11-cis-retinol in the retinal pigment epithelium; due to a membrane-associated alcohol dehydrogenase, 11 cis-retinol is oxidized and converted into 11-cis-retinaldehyde which is the chromophore for rhodopsin and the cone photopigments (By similarity). Required for the survival of cone photoreceptors and correct rod photoreceptor cell morphology (PubMed:25416279). [UniProtKB/Swiss-Prot Function]