

Product datasheet for MR225688

Tpm3 (NM_022314) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tpm3 (NM_022314) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tpm3
Synonyms:	gamma-TM; hTM30nm; hTMnm; Tm5NM; TM30nm; TMnm; Tpm-5; Tpm5; Trop-5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR225688 representing NM_022314 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGGAGGCCATCAAGAAAAAGATGCAGATGCTGAAGTTAGACAAAGAGAATGTTCTGGACCGAGCTG
AGCAAGCTGAAGCTGAGCAGAAACAGGCAGAAGAAAGAAGCAAGCAGCTAGAGGATGAAGTGAAGCAACCAT
GCAGAAGAAGCTGAAAGGGACAGAGGATGAGCTGGACAAGTATTCGGAAGCTTTAAAGGATGCTCAGGAG
AAGCTGGAGCTAGCAGAGAAGAAGGCAGCCGACGCTGAAGCTGAGGTGGCCTCCTTGAACCGCAGGATCC
AGCTGGTTGAAGAGGAGCTGGACCGTGCAGGAGCGCCTTGCCACTGCTTGCAGAAGCTGGAGGAAGC
AGAGAAGCTGCTGATGAGAGTGAGAGAGGTATGAAGGTGATTGAAAATCGGGCTCTAAAAGATGAAGAA
AAGATGGAATCCAGGAAATCCAGCTAAAAGGAAGCAAAGCACATTGCAGAAGAGGCCGATAGGAAGTATG
AAGAGGTGGCTCGTAAGTTGGTATTATTGAAGGAGACTTGGAACGCACGGAGGAACCTGCTGAGCTGGC
AGAGTCCCCTGGCCGAGAGATGGATGAACAGATCAGACTGATGGACCAGAACCTGAAGTGTCTGAGTGTCT
GCTGAAGAAAAGTACTCTCAAAAAGAAGACAAGTATGAAGAAGAAATAAGATTCTTACTGATAAACTCA
AGGAGGCAGAGACCCGTGCTGAGTTTGTGAAAGATCGGTAGCCAAGCTGGAGAAGACCATTGATGACCT
GGAAGATAAACTGAAGTGCACCAAAGAGGAGCATCTCTGTACACAAAGGATGCTGGACCAGACCCCTGCTG
GACCTGAACGAGATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR225688 representing NM_022314
Red=Cloning site Green=Tags(s)


MMEAIKKKMQMLKLDKENVLDRAEQAEAEQKQAEERSKQLEDELATMQKKLKGTEDELDKYSEALKDAQE
 KLELAEKKAADAEAEVASLNRRIQLVEELDRAQERLATALQKLEEAKEAADESERGMKVIENRALKDEE
 KMELQEIQLKEAKHIAEEADRKYEEVARKLVIIIEGDLERTEERAELAESRCREMDEQIRLMDQNLKCLSA
 AEEKYSQKEDKYEEEEIKILTDKLKEAETRAEFAERSVAKLEKTIDDLKDKCTKEEHLCTQRMLDQTL
 DLNEM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9108_d08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_022314


ORF Size: 855 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022314.3 , NP_071709.3
RefSeq Size:	2239 bp
RefSeq ORF:	858 bp
Locus ID:	59069
Cytogenetics:	3 F1
MW:	33.3 kDa
Gene Summary:	Binds to actin filaments in muscle and non-muscle cells. Plays a central role, in association with the troponin complex, in the calcium dependent regulation of vertebrate striated muscle contraction. Smooth muscle contraction is regulated by interaction with caldesmon. In non-muscle cells is implicated in stabilizing cytoskeleton actin filaments.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR225688