

Product datasheet for **MR225711**

Nfatc2 (NM_001037177) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nfatc2 (NM_001037177) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nfatc2
Synonyms:	AI607462; NF-ATc2; NF-ATp; NFAT1; NFAT1-D; Nfatp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR225711 representing NM_001037177
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGGACGTCCCGAGCCGCAGCCCGACCCCGATGGCGGGGACGGCCCCGGCCACGAGCCCGGGGCAGTC
 CCCAAGACGAGCTGGACTTTTCCATCCTCTTCGATTATGACTATCTGAACCCTATCGAAGAAGAACC
 CGCACATAAGGCCATCAGCTCACCTCCGGACTCGCATACCCGGATGATGTCCTGGACTATGGCCTCAAG
 CCATGCAACCCCTTGCCAGTCTCTCTGGCGAGCCCCCTGGCCGGTTCGGAGAGCCGGATAGTATAGGGT
 TCCAGAACTTTCTGAGCCGGTCAAGCCAGCAGGGGCTTCGGGCCGAGCCCTCGGATCGAGATCACTCC
 ATCCCACGAACTGATGCAGGCAGGGGGGCCCTCCGTGGGAGAGACGCCGGCCTGTCCCCGAGCAGCCG
 GCCCTGGCCCTGGCCGGCTGGCCGCCAGCCGAGGTTACACTGCCCGTGCCGGCTACGAGGGCTACC
 GCGAGCCGCTTTGCTTGAGCCCGCTAGCAGCGGCTCCTCTGCCAGTTTCAATTTCTGACACCTTCTCCC
 CTACACCTCGCCCTGCGTCTACCCAATAACGCCGGGCCGACGACCTGTGTCCCGAGTTTCAAACATC
 CCTGCTATTATTCGCCAGAACCTCTCCAATAATGTCACCTCGAACCAGCCTCGCCGAGGACAGCTGCC
 TGGGCCGACACTCGCCCGTGCCCGTCCGGCATCCCGCTCCTCCTCACCCGGTGCCAAGCGGAGGATTC
 GTGCGCAGAGGCTTTGGTTGCTCCTCTGCCCGCAGCCTCACCCAGCGCTCCCGGAGCCCTCGCCACAG
 CCCTCGCTCAGTGGCACTGCAGGACGACAGCATCCCGCTGGGTACCCCCACGCGCCGGCTCTGCTG
 TTCTCATGGATGCCCTCAACACCTGGCCACCGACTCGCCCTGCGGGATCCCTCCAAGATATGGAAGAC
 CAGTCTGACCCGACGCTGTGTCCACCGCTCCGTCCAAGGCTGGCTGGCCGCCACATCTACCTACT
 GTGGAGTTCCTGGGCCATGTGAGCAGGAGGAGAGGAAATTCGCTCCAGAGTCCATCCTGCTGGTAC
 CACCTACTTGGCCAAGCAGTTGGTGGCCGATTCATCTGCAGCATCCCTGTGACTGCATCCCTCCC
 ACCACTCGAGTGGCCACTCTCCAATCAGTCGGGCTCCTATGAGCTACGGATTGAGTCCAACCAAGCCC
 CATCACCGGGCCACTATGAGACGGAGGCGAGCCGTGGCGCTGTCAAAGCCCCAACAGGAGGACACCCTG
 TGGTGCAGCTCCACGGCTACATGGAGAACAAGCCTCTGGGGCTTCAGATCTTCATTGGGACAGCAGATGA
 GAGGATCCTTAAGCCGCACGCTTCTACCAAGTACACAGGATCACTGGGAAAACGGTACCACCACGAGC
 TATGAGAAGATCGTAGGCAACACCAAGTCTGGAGATCCCCCTGGAGCCAAAGAACAACATGAGAGCCA
 CCATCGACTGTGCAGGCATCCTGAAGCTCCGAAACGCTGACATCGAGCTGCGGAAGGGCAGACGGACAT
 CGGCAGGAAGAACACGCGTGTGCGCCTGGTGTCCGCGTGCACGTCCAGAGCCAGTGGGCGCATCGTC
 TCCTGCAGGCTGCGTCAACCCATCGAGTGTCTCAGCGCTCTGCCACGAGCTGCCATGGTGGAGA
 GACAAGACATGGACAGCTGCCTGGTCTACGGGGCCAGCAGATGATCCTCACGGGCCAGAACTTCACAGC
 GGAGTCCAAGGTTGTGTTTCATGGAGAAGACTACAGGACCTGCAGGGACCTGTGAGACTCGGCCTTTGCC
 ATCTCTGATCTCAGCTGACCGTCTCTCCCCCTGGCTTTCCCGGCTGCAGAGAAACCTCCTGGCTCTG
 TCTTCAGATGCTCCGTGCTCCTCCTGCCCCAGGATCTTCGCTTGTGCTGCTGGCTCTG

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR225711 representing NM_001037177
Red=Cloning site Green=Tags(s)


```
MDVPEPQPPDGGDGPHEPGGSPQDELDF SILFDYDYLNPIDEEPIAHKAISSPSGLAYPDDVLDYGLK
PCNPLASLSGEPGRFGEPDSIGFQNFLSPVKPAGASGSPRIETPSHELMQAGGALGRDAGLSPEQP
ALALAGVAASPRFTLPVPGYEGYREPLCLSPASSGSSASFISDTFSPYTPSPCVSPNAGPDDLCPQFQNI
PAHYSPRTSPIMSPRTSLAEDSCLGRHSPVPRPASRSSSPGAKRRHSACAEALVAPLPAASQRSRSPSPQ
PSPHVALQDDSI PAGYPPTAGSAVLM DALNTLATDSPCGIPSKIWKTS PDPTPVSTAPSKAGLARHIYPT
VEFLGPCEQEERRNSAPESILLVPPTWPKQLVPAIPICISIPVTASLPLEWPLSNQSGSYELRIEVQPKP
HHRAHYETEGSRGAVKAPTGGHPVVQLHG YMENKPLGLQIFIGTADERILKPHAFYQVHRITGKVTTTTS
YEKIVGNTKVLEIPLPKNNMRATIDCAGILKLRNADIELRKGETDIGRKNTRVRLVFRVHVPEPSGRIV
SLQAASNPIECSQRSAHELPMVERQDMDSCLVYGGQQMILTGQNF TAESKVVFMKTTGPAGTCETRPLP
ISLISADRLSPWLSRLQRNPPGSVFRCSVLLPAPGSSLVLLAL
```

SGPTRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-RsrII

Cloning Scheme:



ACCN: NM_001037177

ORF Size: 2019 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001037177.2](#), [NP_001032254.1](#)

RefSeq Size: 3362 bp

RefSeq ORF: 2022 bp


Locus ID: 18019

Cytogenetics: 2 88.91 cM

MW: 72.4 kDa

Gene Summary: Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2, IL-3, IL-4, TNF-alpha or GM-CSF. Promotes invasive migration through the activation of GPC6 expression and WNT5A signaling pathway.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR225711