

Product datasheet for **MR225738L4V**

Pard3 (NM_001013580) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pard3 (NM_001013580) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pard3
Synonyms:	AA960621; AI256638; Asip; D8Ertd580e; Par-3; Par3; Pard-3; Pard3a; Phip
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001013580
ORF Size:	1818 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225738).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001013580.3
RefSeq Size:	3629 bp
RefSeq ORF:	1821 bp
Locus ID:	93742
UniProt ID:	Q99NH2
Cytogenetics:	8 74.66 cM



[View online »](#)

Gene Summary:

Adapter protein involved in asymmetrical cell division and cell polarization processes (By similarity). Seems to play a central role in the formation of epithelial tight junctions (By similarity). Targets the phosphatase PTEN to cell junctions (By similarity). Association with PARD6B may prevent the interaction of PARD3 with F11R/JAM1, thereby preventing tight junction assembly (PubMed:11839275). The PARD6-PARD3 complex links GTP-bound Rho small GTPases to atypical protein kinase C proteins (By similarity). Required for establishment of neuronal polarity and normal axon formation in cultured hippocampal neurons (By similarity). Involved in Schwann cell peripheral myelination (PubMed:21949390). [UniProtKB/Swiss-Prot Function]