

Product datasheet for **MR225749**

Galc (NM_008079) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Galc (NM_008079) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Galc
Synonyms:	2310068B06Rik; A930008M05Rik; AW212969; AW413532; Gac; Gacy; twi; twitcher
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR225749 representing NM_008079
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTAACAGCCAACCTAAGGCTTCCCAGCAACGCCAAGCAAAAGTCATGACCGCCGCCGCGGGCTCGG
 CGAGCCGTGTTGCGGTGCCCTTATTGTTGTGTGCGCTGCTAGTGCCCGGTGGCGCCTACGTGCTAGACGA
 CTTCCGACGGGCTGGGCCGGGAGTTTCGATGGCATCGGCGCAGTCAGCGCGCGGGCAACCTCCCGGCTT
 CTAGTAAATTACCCAGAGCCCTATCGTTCTGAAATACTGGATTATCTATTTAAGCCAAACTTTGGTGCTT
 CTTTGCATATTCTAAAAGTGGAGATAGGTGGTATGGACAAACAACAGATGGCACTGAACCTCGCACAT
 GCACTATGAGCTAGATGAGAATTATTTCCGAGGCTATGAGTGGTGGCTAATGAAGGAAGCTAAGAAGCGG
 AACCCAGATATTACTTATGGGTTGCCATGGTCATTCCCTGGATGGCTGGGCAAAGGCTTCAGCTGGC
 CTTATGTAATCTTCAGCTCACTGCCTACTATGTTGTAAGATGGATTTTGGGTGCCAAGCATTATCATGA
 CCTGGACATTGATTATATTGGGATTTGGAATGAAAGCCGTTTGACGCCAACTATATAAAGGAATTAAGA
 AAAATGCTGGATTATCAAGGTCTCCAGCGAGTGAGAATCATAGCGAGCGATAATCTCTGGGAGCCGATTT
 CCTCTTCCTTGCTGCTGGACCAGGAACCTGGAAGGTGGTTGATGTTATAGGAGCTCACTATCCTGGAAC
 CTACACAGTGTGGAATGCAAAGATGTCAGGGAAGAAGCTGTGGTCATCTGAAGATTTTAGCACTATCAAC
 AGTAATGTTGGCGCAGGCTGTTGGAGTCGCATATTGAACCAGAATTACATCAATGGCAATATGACCTCCA
 CAATTGCTTGGAAATTTGGTGGCTAGTACTATGAAGAGTTGCCTATGGACGAAGTGGGTTGATGACGGC
 CCAGGAGCCATGGAGTGGACATTATGTCGTAGCATCTCCTATCTGGGTATCAGCTCACACAACCCAGTTT
 ACTCAACCAGGCTGGTATTACCTGAAGACAGTGGGCCATCTAGAGAAAGGAGGAAGTTATGTTGCTCTGA
 CTGATGGCTTGGGAAACCTCACCATCATCATTGAAACTATGAGTCATCAACATTCTATGTGTATACGGCC
 ATATCTGCCATATTATAATGTCTCACATCAGCTGGCTACCTTCACTCTAAAGGGATCTTTGAGAGAATA
 CAAGAGCTGCAGGTGTGGTATACCAAGCTGGGAACTCCGCAGCAGAGACTTCATTTTAAACAACCTGGATA
 CTCTATGGCTCCTTGACGGCAGCGGAGTTTCACTGGAACCTGGAGGAGGACGAGATATTTACGCTCAC
 CACTCTGACCACAGGTCGCAAAGGCAGCTACCCTCCCCCTCCAGCTCGAAGCCCTTCCCAACCAACTAC
 AAGGATGATTTCAATGTTGAGTACCCACTTTTTAGTGAAGCTCCGAATTTTGTGATCAGACTGGCGTGT
 TTGAGTACTACATGAATAATGAAGACCGTGAGCACCGCTTACGCTGCGCCAAGTTCTCAACCAACGACC
 TATTACCTGGGCTGCAGACGCTCCAGCACAATCAGTGTATAGGCGATCACCCTGGACCAACATGACT
 GTGCAGTGGATGTTTATATAGAGACCCCTCGGAGCGGCGGCGTATTCTTGTGGAAGAGTGAATAAAG
 GTGGAATCTTGATTAGAAGTGTACAGGAGTCTTTTTCTGGATTTTGGCAATGGCTCTTACAGAGTCAC
 AGCTGACCTGGGTGGATGGATAACGTATGCTTCAGGACATGCGGACGTTACAGCTAAACGATGGTACACA
 CTCACCTTAGGCATTAAGGGTATTTCCGCTTTGGTATGTTGAACGGCACGATCCTATGGAAAAACGTCC
 GTGTGAAATACCCAGGCGATGGCTGGGCTGCGATCGGAACACACACCTTCGAGTTTGCAGTTTGACAA
 CTTTCGCGTGGAAAGCTGCTCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_008079.4 , NP_032105.2
RefSeq Size:	3723 bp
RefSeq ORF:	2055 bp
Locus ID:	14420
UniProt ID:	P54818
Cytogenetics:	12 49.83 cM
MW:	77.7 kDa
Gene Summary:	This gene encodes galactosylceramidase, the lysosomal hydrolase involved in the catabolism of galactosylceramide. Mutations in this gene result in slow growth, tremors and hind leg weakness, collectively termed as the 'twitcher' phenotype. In humans, deficiency of this gene product causes a lysosomal storage disorder known as Krabbe disease. [provided by RefSeq, Dec 2014]

Product images:



Circular map for MR225749