

## Product datasheet for **MR225763L4V**

### Hipk3 (NM\_010434) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hipk3 (NM_010434) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hipk3
Synonyms:	DYRK6; FIST3; mir-1902; PKY
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010434
ORF Size:	3576 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225763).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010434.2</a> , <a href="#">NP_034564.2</a>
RefSeq Size:	7497 bp
RefSeq ORF:	3579 bp
Locus ID:	15259
UniProt ID:	<a href="#">Q9ERH7</a>
Cytogenetics:	2 E2



[View online »](#)

**Gene Summary:**

Serine/threonine-protein kinase involved in transcription regulation, apoptosis and steroidogenic gene expression. Phosphorylates JUN and RUNX2. Seems to negatively regulate apoptosis by promoting FADD phosphorylation. Enhances androgen receptor-mediated transcription. May act as a transcriptional corepressor for NK homeodomain transcription factors. The phosphorylation of NR5A1 activates SF1 leading to increased steroidogenic gene expression upon cAMP signaling pathway stimulation. In osteoblasts, supports transcription activation: phosphorylates RUNX2 that synergizes with SPEN/MINT to enhance FGFR2-mediated activation of the osteocalcin FGF-responsive element (OCFRE).[UniProtKB/Swiss-Prot Function]