

## Product datasheet for **MR225806L3V**

### Meis1 (NM\_001193271) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Meis1 (NM_001193271) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Meis1
Synonyms:	C530044H18Rik; Evi8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001193271
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225806).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001193271.1</a> , <a href="#">NP_001180200.1</a>
RefSeq Size:	3441 bp
RefSeq ORF:	1173 bp
Locus ID:	17268
UniProt ID:	<a href="#">Q60954</a>
Cytogenetics:	11 11.11 cM



[View online »](#)

**Gene Summary:**

Acts as a transcriptional regulator of PAX6. Also acts as a transcriptional activator of PF4 in complex with PBX1 or PBX2. Required for hematopoiesis, megakaryocyte lineage development and vascular patterning. May function as a cofactor for HOXA7 and HOXA9 in the induction of myeloid leukemias.[UniProtKB/Swiss-Prot Function]