

Product datasheet for **MR225808**

Gpr155 (NM_001080707) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpr155 (NM_001080707) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gpr155
Synonyms:	1110017O10Rik; AW319992; DEPDC3; F730029F15Rik; PGR22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR225808 representing NM_001080707
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGATTCTTATTCTCAGCAAAGAACTCAACGCTTGCTGGTGATATGAATGCGACTTGGCTGCCTCTC
 ATGGATTTAATGCCACCGCGATCCACCTTCGATGTCGATAACCCGGCTTCCCGGCCTTACTGGAGTG
 CTTTGGCATAGTCCTTTGTGGCTACATAGCAGGAAGGGCCAATATCATAACCTCAACACAGGCCAAAGGA
 CTGGGGAAATTCGTCCTCAGATTTGCCCTTCCAGCCTTATTATTCAAAAACATGGTTGTGCTTAATTTTT
 CCAACGTGGACTGGGCTTTTCTACAGCGTCTTAATTGGCAAAGCCTCCGTGTTTTTCATTGTGTGTGT
 ATTAACCTTATTGGTGGCCAGTCTGAGAGCCGGTTTAGCAAAGCTGGACTATCCCTATTTTTGCTACA
 CAGAGTAATGACTTTGCTTTGGGATACCCTATAGTTGAAGCTCTGTACCAATCTACATATCCAGAGTATC
 TGCAGTACATTTATTTGGTGGCACCAATATCTCTTATGATGCTAAACCCCATAGGGTTTATCTTCTGTGA
 AATCCAGAAGTCAAAGACACTCAAATGCTTCTCAAATAAAGCAAAAATTGTGGGACTTGGATTTCTA
 CGCGTGTACAGAACCCGATAGTATTTATGGTCTTTGTTGGCATTGCCTTCAATTTTATTCTTGATAAAA
 AGATCCAGTATATATGGAAAACCTTCTGGATGGGCTTCAAATTCCTTCTCTGGATCAGCCCTGTTTTA
 TCTTGGTCTGACCATGGTAGGAAAAATCAGGCGACTGAAGAAGTCGGCGTTTGTGGTTCTTACTCTTCTC
 ATCACAGCCAACTCCTGGTGTGCCCTCCTGTGCAGAGAAATGGTGGAACTCCTGGACAAGGGGGACA
 GCGTAGTGAACCATACAAGCCTGTCCAATATGCGTTTCTATATGGTGTCTTCCCTGTGGCTCCAGGAGT
 GGCTATCTTTGCAACACAGTTCAACATGGAAGTGGAGATTATAACCTCAGGCATGGTAATAAGTACATTT
 GTGTCTGCTCCAATCATGTATGTTTCTGCCTGGTTATTAACCTTCTACCATGGATGCCAAGCCACTGG
 CGTATGCCAATCAGAATGTTAGCTTTGATATAAGTATTATCAGCCTGGTCTCCTTGATCTGGTCCCTGTC
 TATTTCTCTTGAGTAAGAAAATAAGCAGCTCCACATATGCTCACAGCGAACTTACTCATCGCTCAG
 ACTATTGTTTGTGCTGGAATGATGATATGGAATTTTGTAAAGAAAAAATTTTGTGGACAAAATTTGG
 TGTTTTGTTTTATTGTACAGCTCCCTGTACAGCACCTATCTGTGGACAGGCCTTCTAGCAGTTTCTTTGTT
 TCTTTAAAGAAGCGGAGAGCGTTCAGCTCCCGGTGGGAATCATCATCATCTCTGGCTGGGGGATCCCT
 GCTCTCCTTGTGGCGTCTTCTGATAACTGGAACACACACGGAGACAGCATTGACTCAGCCTTCTTTT
 ATGAAAAGAGCAGATGATCACACGGCAGTACCCTGTTCTGCAGATTCTCATAGCTGGGGTGTCACT
 CATGTGCATGAACCGCACCAAGCTGGTCACTACGAAGTTTCGGCCAGTCTCAGAACCACAAACCA
 GTAGAGCCTGGAAGCACTGCGTTTGGAGAGAACCAGCACCAACAAATGAACCAGAACTTCCCAAGTT
 CTATCCAGAAAAGTTGCTGTTCTTGTCTCTTGGGAAATGGTGAATTACGCTGCCCATCCATAGAGCC
 AGTAACGAACTCAAGTGCTAGTGGGCTATGCCTTCTTCTTTGAGAAAACCGACCACTGTGTGAGCCGT
 TGTGACTCCCAGAGCTGCATCCTCGCCAGGAAGAAGAGCAGTACCTGCAGAGCGGAGACCCACAGCTGA
 CCCGGCACGTCTGTTGTGCTTGTCTTATCATCGGCCTCTTTGCTAATCTCTCCAGCTGCCTGTGGTG
 GTTGTTCACACGAGACGGGACGGCTCTATGTTGAATTACAGTTTTTCTGTGCCGTGTTAACTTTGGC
 CAGTGGGTGCTGCGTAT

AG**GCGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR225808 representing NM_001080707
Red=Cloning site Green=Tags(s)

```
MDSYFSAKNSTLAGDMNATWPASHGFNATGDPPSMSITRLFPALLECFGIVLCGYIAGRANIITSTQAKG
LGNFVSRFALPALLFKNMVVLFNSVDWAFLYSVLIGKASVFFIVCVLTLLVASPESRFKAGLFPFIAT
QSNDFALGYPIVEALYQSTYPEYLQYIYLVAPISLMLNPIGFIFCEIQKSKDTQNASQNKAKIVLGLFL
RVLQNPVIVFMVFGIAFNFIIDKKIPVYMFENFLDGLANSFSGSALFYGLTMTVGKIRRLKKSFAFVLTLL
ITAKLLVLPLLCREMVELLDKGDVSVNHTSLSNYAFLYGVFPVAPGVAIFATQFNMEVEIITSGMVISTF
VSAPIMYVSAWLLTFPTMDAKPLAYAIQNVSFDISIISLVSLIWSLSILLSSKKYKQLPHMLTANLLIAQ
TIVCAGMMIWNFVKEKNFVQILVFLVLYSSLYSTYLWTGLLAVSLFLLKRESVQLPVGIIISGWGIP
ALLVGVLLITGKHNGDSIDSAFFYGKEQMITTAVTLFCSIL IAGVSLMCMNRTTQAGHYEGFGSQNHKP
VEPGSTAFEENPAPTNEPELFPSSIPETSCCSCSLGNGELRCPSEIPTVNTSSASGPMPSSEFKTDHCVSR
CDSQSCILAQEEEQYLQSGDPQLTRHVLLCLLLIIGLFANLSSCLWWLFNHETGRLYVELQFFCAVFNFG
QVGAAY
```

SGPTRRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001080707

ORF Size: 2118 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080707.1](#), [NP_001074176.1](#)

RefSeq Size: 2625 bp


RefSeq ORF: 2121 bp

Locus ID: 68526

Cytogenetics: 2 C3

MW: 96.3 kDa

Product images:



Circular map for MR225808