

## Product datasheet for **MR225840**

### **Pard3 (NM\_001013581) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Pard3 (NM_001013581) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pard3
Synonyms:	AA960621; AI256638; Asip; D8Ertd580e; Par-3; Par3; Pard-3; Pard3a; Phip
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR225840 representing NM\_001013581  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**


ATGCCTCTTCATGTCCGCCGGAGCAGCGACCCAGCTTTAACTGGCCTTTCCACTTCTGTCTAGTGATAACA  
 ACTTTTCTCAGAGGAGCCCTCCAGGAAAAACCCACCCGCTGGTCCACGACAGCTGGCTTTCTCAAGCA  
 GAACACCGCTGGAAGTCCCAAAACCTGCACAGGAAGAAAGATGAAAACCTACAGAAGCCTTCCACGGGAT  
 CCCAGTAGCTGGTCCAACCAAGTCCAGCGAGACAACGCCCGCTCCTCCCTGAGCGCCAGCCACCCAAATGG  
 TAGACCGGTGGCTGGAGAAGCAAGAACAGGATGAGGAAGGCACAGAAGAAGACAGCAGCCGAGTGGAGCC  
 GGTTGGACATGCTGATACCGGATTGGAGAACATGCCCAACTTTTCCCTCGATGATATGGTAAAGCTCGTA  
 CAAGTCCCAACGATGGAGGGCCCTGGGAATCCATGTAGTGCCTTTTCAGTGTCTGAGCGCCGAGAACAT  
 TGGGGTTGTTAGTGAAGCGGTTGGAGAAAGCGGTAAGGCTGAGCAAGAAAACCTTTTCCATGAGAATGA  
 CTGCATTGTGAGGATTAACGATGGAGATCTTCGAAACAGAAGATTTGAGCAAGCACAACATATGTTCCGC  
 CAAGCTATGCGTGCAGCTGATTTGGTCCATGTGGTCCCTGCAGCAAACAAGGAGCAATATGAACAAC  
 TGTCCTCAACGCGAGAAGAACAACACTACTCCCAAGGCCGCTTCAGCCCTGACAGCCACTGTGTGGCCAAACAG  
 GAGTGTGGCCAACAATGCCCTCAAGCATTGCCAGAGACCCAGACTGAGTCAAGCCACCCGAGCAGCTG  
 GATGCTCACCCCGACTACCTCATAGTGCTCACGCCTCAACCAAACCCCGCAGCCCGGCCCTTGCTC  
 CACCCAGTGTGCTTAGTACCAACGTAGGCAGTGTGTACAACACGAAGAAAGTAGGCAAGAGGCTCAACAT  
 CCAGCTTAAGAAAGGTACAGAAGGACTGGGATTCAGCATCACCTCCCGGACGTCACCATAGGTGGCTCA  
 GCTCCCATTTATGTCAAGAATATCCTTCTCGAGGGGCTGCCATTCAGGATGGCAGACTCAAGGCAGGAG  
 ACCGGCTAATAGAGGTCAATGGAGTAGATTTAGCAGGCAAAATCCAGGAGGAAGTTGTTTCCCTGTTGAG  
 AAGCACCAAGATGGAGGGGACTGTGAGCCTTCTGGTCTTTCGTCAGGAAGAGGCTTCCACCCAAGGGAA  
 ATGAATGCTGAACCAAGCCAGATGCAGACTCCTCAAAAGAAACGAAAGCTGAAGATGAGGACGTTGTTCTCA  
 CACCCGATGGTACCAGGGAGTTTCTGACTTTTGAAGTTCCACTGAATGACTCAGGATCTGCAGGGCTTGG  
 TGTCAAGTGTCAAGGGGAACCGTTTCAAAGAGAACCACGCAGATTTGGGGATCTTCGTTAAATCCATTATC  
 AATGGTGGAGCTGCATCTAAGATGGAAGGCTGAGGGTAAATGACCAGCTGATAGCTGTGAATGGAGAAT  
 CTCTACTGGGCAAAGCCAACCAGGAAGCCATGGAGACTCTACGGAGGTCCATGTCCACCGAGGGCAACAA  
 GCGTGGCATGATCCAGCTCATTGTGGCGAGGCGGATCAGCAGATGTAACGAGCTGCGGTCTCCTGGGAGC  
 CCTGCTGCACCTGAGCTGCCATCGAGACAGAAGTGGATGACCGAGAACGCAGGATCTCACACTCCCTCT  
 ACAGTGGGATCGAGGGGCTGGATGAGTCGCCCACCAGGAACGCAGCACTGAGCAGGATAATGGGTAATG  
 CCAGCTGTCCCCACGGTGAACATGCCTCACGATGACACTGTCTGATTGAAGATGACAGGCTGCCTGTG  
 CTCCCTCCTCACCTCTCTGACCAGTCTCCTCCAGCTCCCATGATGACGTTGGATTGATAATGACAGAAG  
 CAGGCAGTGGGCCAAAGCTACCATCAGTGACTCAGCCGACTGCTCATTGAGTCCAGATGTTGATCCGGT  
 TCTCGTTTTCAACGGGAAGGATTTGGACGCCAGAGTATGTCAGAAAAACGCACAAGCAATTTTCAGAT  
 GCCAGTCAATTGGATTTTCGTTAAAACACGAAAATCAAAAAGCATGGATTTAGGTATAGCTGACGAGACTA  
 AACTCAATACAGTGGATGACCAGAGAGCAGGCTCTCCAGTAGAGATGTGGGACCCTCCTTAGGTCTGAA  
 GAAATCTAGCTCCTTAGAAAGTCTGACAGCCGCTGTCGCGAGGTGACGCTGAACGGGAACATTCCTTTT  
 CACCGCCCTCGGCCACGAATCATCCGGGAAGGGGCTGCAACGAAAGCTTCAGAGCCGCCATTGACAAGT  
 CCTACGACAAGCCATGTTGATGATGACGACGAAGGCATGGAGACCTTGAAGAAGACACAGAGGAAAG  
 CTCGAGGTGAGGGAGGGAGTCCGTGTCCAGTCCAGCGATCAGCCTTCTATTCTCTGGAGAGGCAATG  
 AATGGAGACCCAGAGAAAAGGGACAAGACAGAGAGGAAAAAGGACAAAGCCGAAAGGATAAGAAGAAAG  
 ACCGAGAGAAGGAGAAGGATAAACTGAAAGCCAAGAAGGGGATGCTGAAAGGCTTGGGGGACATGTTTCAG  
 CCTGGCCAAACTGAAGCCGGAGAAGAGA

**ACCGGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



<b>ORF Size:</b>	2688 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001013581.3</a>
<b>RefSeq Size:</b>	3233 bp
<b>RefSeq ORF:</b>	2691 bp
<b>Locus ID:</b>	93742
<b>UniProt ID:</b>	<a href="#">Q99NH2</a>
<b>Cytogenetics:</b>	8 74.66 cM
<b>MW:</b>	99.3 kDa
<b>Gene Summary:</b>	Adapter protein involved in asymmetrical cell division and cell polarization processes (By similarity). Seems to play a central role in the formation of epithelial tight junctions (By similarity). Targets the phosphatase PTEN to cell junctions (By similarity). Association with PARD6B may prevent the interaction of PARD3 with F11R/JAM1, thereby preventing tight junction assembly (PubMed:11839275). The PARD6-PARD3 complex links GTP-bound Rho small GTPases to atypical protein kinase C proteins (By similarity). Required for establishment of neuronal polarity and normal axon formation in cultured hippocampal neurons (By similarity). Involved in Schwann cell peripheral myelination (PubMed:21949390). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR225840