

Product datasheet for **MR225865**

Keap1 (NM_016679) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Keap1 (NM_016679) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Keap1
Synonyms:	INRF2; mKIAA0132
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR225865 representing NM_016679
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCAGCCGAACCCAAGCTTAGCGGGGCTCCCGCAGCAGCCAGTTCCTGCCCTGTGGTCAAAGTGCC
 CCGAGGGGGCCGGGACGCAGTGATGTATGCCCTCCACGGAGTGCAAGGCAGAGGTGACGCCCTCGCAGGA
 CGGTAACCGAACCTTCAGCTACACACTAGAGGATCACACCAAGCAGGCTTTTGGCGTCATGAACGAGCTT
 CGCCTGAGCCAGCAACTCTGTGACGTGACCCTGCAGGTCAAATATGAGGACATCCAGCTGCCAATTCA
 TGGCTCACAAAGTGGTGTGGCCTCTCCAGCCAGTCTTTAAAGCCATGTTACCAACGGGCTTCGGGA
 GCAGGGCATGGAGGTGGTGTCCATCGAAGGCATCCACCCTAAGGTATGGAAAGGCTTATTGAGTTCCGC
 TACACGGCCTCCATCTCCGTGGCGGAGAAGTGTCTGCACGTGATGAACGGGGCGGTATGTACCAGA
 TTGACAGCGTGGTTTCAGCCTGCAGCGACTTCTCGTGCAGCAGCTGGACCCAGCAACGCCATTGGCAT
 CGCCAACTTCGCGGAGCAGATCGCTGCACTGAACTGCACCAGCTGCCCGGAGTATATCTACATGCAC
 TTCGGGGAGGTGGCAAGCAGGAGGAGTCTTCAACCTGTCACTGCCAGCTGGCCACGCTCATCAGCC
 GGGATGATCTGAACGTACGCTGCGAGTCCGAGGTGTTCCACGCGTGCATCGACTGGGTCAAATACGACTG
 CCCGCAGCGGCGCTTCTACGTGCAGGCACTGCTGCGGGCGTGCCTGCCATGCGCTCACGCCGCGCTTC
 CTGCAGACACAGCTGCAGAAAGTGTGAGATCCTGCAGGCCGACGCGCGTGCAGGACTACCTGGTGCAGA
 TATTCCAGGAGCTCACGTGCACAAGCCACGCAGGCACTGCCCTGCCGCGGCCCAAAGTGGGCCGCT
 CATCTACACAGCGGGCGGTTACTTCCGACAGTCGCTCAGCTACCTGGAGGCTACAACCCGAGCAATGGC
 TCCTGGCTGCGCTGGCCGATCTACAGGTGCCGCGCAGTGGGCTGGCAGGCTGCGTGGTGGTGGGCTGC
 TATACGCTGTGGCGGCGCAACAACACTCCGGATGGCAACACTGACTCCAGCGCCCTGGACTCTACAA
 CCCCATGACCAACCAGTGGTCCGCTGTGCCTCTATGAGCGTGCCACGCAACCGCATCGGGGTGGGGGTC
 ATAGATGGCCACATCTACGAGTCGGGGGTTCCACGGCTGCATCCACCACAGCAGCGTGGAGAGATATG
 AGCCAGAGCGGGACGAGTGGCATCTAGTCGCGCAATGTTGACACGGAGGATTGGCGTGGGCGTGGCAGT
 GCTCAACCGCTTGTGTATGCAGTGGGGGCTTTGACGGGACTAACCGGCTTAACTCCGAGAAATGTTAC
 TATCCAGAGAGGAATGAGTGGCGGATGATCACACCGATGAATACCATCCGAGCGGGGCGGGGCTGCGG
 TGCTGCACAACGTATCTATGCAGCAGGGGCTACGATGGCAGGACCAGTTGAACAGTGTGGAGCGCTA
 CGAGTGGAGACAGAGACCTGGACTTTCGTAGCCCCATGAGGCATCACCGTAGTGGCTGGGGATTACT
 GTGCACCAGGGCAAGATCTACGTCTCGGAGGCTATGATGGCCACACTTTCTGGACAGTGTGGAATGCT
 ATGACCCGGACAGTGATACCTGGAGTGAGGTGACCCGATGACATCTGGCCGACGCGGGTGGGTGTGGC
 CGTCACCATGGAACCCTGTCGGAAGCAAATTGATCAACAAAACGTACCTGC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR225865 representing NM_016679
 Red=Cloning site Green=Tags(s)

MQPEPKLSGAPRSSQFLPLWSKCEPAGDAVMYASTECKAEVTPSQDGNRTFSYTLIEDHTKQAFGVMNEL
 RLSQQLCDVTLQVKYEDIPAAQFMAHKVVLASSPVFKAMFTNGLREQGMEVVSIEGIHPKVMERLIEFA
 YTASISVGEKCVLHVMNGAVMYQIDSVVRACSDFLVQQLDPSNAIGIANFAEQIGCTELHQRAREIYMH
 FGEVAKQEEFFNL SHCQLATLISRDDLNVRCSEVVFHACIDWVKYDCPQRRFYVQALLRAVRCHALTPRF
 LQTQLQKCEILQADARCKDYLVQIFQELTLHKPTQAVPCRAPKVGRLIYTAGGYFRQSLSYLEAYNPSNG
 SWLRLADLQVPRSLAGCVVGGLLYAVGGRNNSPDGNTDSSALDCYNPMTNQWSPCAMSVPNRIRIGVGV
 IDGHIYAVGGSHGCIHSSVERYEPERDEWHLVAPMLTRRIGVGVAVLNRLLYAVGGFDGNTNRLNSAECY
 YPERNEWRMITPMNIRSGAGVCLHNCIYAAGGYDGGDQLNSVERYDVETETWTFVAPMRHRSALGIT
 VHQGIYVVLGGYDGHFLDSVECYDPPSDTWSEVTRMTSGRSGVGVAVTMEPCRKQIDQQNCTC

SGP**TRRRLEQKLI**SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_016679

ORF Size: 1872 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016679.1](#), [NM_016679.2](#), [NM_016679.3](#), [NM_016679.4](#), [NP_057888.1](#)
RefSeq Size: 4470 bp
RefSeq ORF: 1875 bp
Locus ID: 50868
UniProt ID: [Q9Z2X8](#)
Cytogenetics: 9 A3
MW: 70 kDa

Gene Summary: Acts as a substrate adapter protein for the E3 ubiquitin ligase complex formed by CUL3 and RBX1 and targets NFE2L2/NRF2 for ubiquitination and degradation by the proteasome, thus resulting in the suppression of its transcriptional activity and the repression of antioxidant response element-mediated detoxifying enzyme gene expression. Retains NFE2L2/NRF2 and may also retain BPTF in the cytosol. Targets PGAM5 for ubiquitination and degradation by the proteasome (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR225865