

Product datasheet for **MR225875L3V**

Fut7 (NM_001177367) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fut7 (NM_001177367) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fut7
Synonyms:	AI853193; FTVII; Fuc-TVII; FucT-VII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177367
ORF Size:	1029 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225875).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001177367.1 , NP_001170838.1
RefSeq Size:	1849 bp
RefSeq ORF:	1029 bp
Locus ID:	14347
UniProt ID:	Q11131
Cytogenetics:	2 A3
Gene Summary:	Catalyzes alpha-1,3 glycosidic linkages involved in the expression of sialyl Lewis X antigens. [UniProtKB/Swiss-Prot Function]



[View online »](#)