

## Product datasheet for **MR225913L3V**

### **Fgf17 (NM\_008004) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Fgf17 (NM_008004) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fgf17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008004
ORF Size:	651 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225913).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008004.4</a>
RefSeq Size:	1687 bp
RefSeq ORF:	651 bp
Locus ID:	14171
UniProt ID:	<a href="#">P63075</a>
Cytogenetics:	14 36.4 cM
Gene Summary:	Plays an important role in the regulation of embryonic development and as signaling molecule in the induction and patterning of the embryonic brain. Required for normal brain development.[UniProtKB/Swiss-Prot Function]



[View online »](#)