

Product datasheet for **MR225919L3V**

Rps6 (NM_009096) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rps6 (NM_009096) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rps6
Synonyms:	S6R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009096
ORF Size:	747 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225919).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009096.3 , NP_033122.1
RefSeq Size:	1394 bp
RefSeq ORF:	750 bp
Locus ID:	20104
UniProt ID:	P62754
Cytogenetics:	4 C4
Gene Summary:	May play an important role in controlling cell growth and proliferation through the selective translation of particular classes of mRNA.[UniProtKB/Swiss-Prot Function]



[View online »](#)