

Product datasheet for **MR225926L4V**

Rbp4 (NM_011255) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rbp4 (NM_011255) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rbp4
Synonyms:	Rbp-4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011255
ORF Size:	603 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225926).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011255.3 , NP_035385.1
RefSeq Size:	930 bp
RefSeq ORF:	606 bp
Locus ID:	19662
UniProt ID:	Q00724
Cytogenetics:	19 32.75 cM
Gene Summary:	Retinol-binding protein that mediates retinol transport in blood plasma. Delivers retinol from the liver stores to the peripheral tissues. Transfers the bound all-trans retinol to STRA6, that then facilitates retinol transport across the cell membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)