

Product datasheet for **MR225939L3V**

Gfy (NM_001195255) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gfy (NM_001195255) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gfy
Synonyms:	Gm581; Gm2515; Goofy
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001195255
ORF Size:	1512 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225939).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001195255.1 , NP_001182184.1
RefSeq Size:	1877 bp
RefSeq ORF:	1515 bp
Locus ID:	100039953
UniProt ID:	J3KML8
Cytogenetics:	7
Gene Summary:	Required for proper function of the olfactory system. May be involved in establishing the acuity of olfactory sensory signaling.[UniProtKB/Swiss-Prot Function]



[View online »](#)