

Product datasheet for **MR225956A1V**

Mouse 2310002J21Rik (NM_001195348) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse 2310002J21Rik (NM_001195348) AAV Particle
Tag: Myc-DDK
Symbol: 2310002J21Rik
Synonyms: 2310002J21Rik; 2610002B11Rik; 4930560M04Rik; Ncrna00081
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR225956 representing NM_001195348
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCAGAAGTGAAGTCGATGTTCCGAGAAGTCCTTCCAAAACAAGGGCAGTTGTCTGTGGAAGATGTAA
 CCACAATGGTGTGTGTAACCCAAGCTCTTGCCTTAAATCTCTGACGCTGGAAAACTGGAGAAAAT
 GCAGCAAGCAGCACAGGATACAGTCCGTCAACAAGAAATGGCAGAAAAGGCAACGAGCAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR225956 representing NM_001195348
 Red=Cloning site Green=Tags(s)

MAEVKSMFREVLPKQQLSVEDVTTMVLCKPKLLPLKSLTLEKLEKMQAAQDTRVQQEMA EKATSK

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001195348
ORF Size: 201 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: NM_001195348.1, NP_001182277.1

RefSeq Size: 2248 bp

RefSeq ORF: 204 bp

Locus ID: 66360

Cytogenetics: 19 D2

MW: 7.6 kDa