

Product datasheet for **MR226000L4V**

Skor2 (NM_001109743) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Skor2 (NM_001109743) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Skor2
Synonyms:	Corl2; EG664805; Fussel18; Gm7348
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001109743
ORF Size:	3024 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226000).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001109743.1 , NP_001103213.1
RefSeq Size:	3141 bp
RefSeq ORF:	3027 bp
Locus ID:	664805
UniProt ID:	A7M7C7
Cytogenetics:	18 E3
Gene Summary:	Acts as a TGF-beta antagonist in the nervous system (By similarity). Exhibits transcriptional repressor activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)