

Product datasheet for **MR226002A1V**

Mouse *Triqk* (NM_001171801) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse *Triqk* (NM_001171801) AAV Particle
Tag: Myc-DDK
Symbol: *Triqk*
Synonyms: C130086A10; Gm11818; OTTMUSG00000004461
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR226002 representing NM_001171801
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGTAGAAAAGACTCTTCCAATACAAAACCTCCTGTTGATCAGTACAGAAAGCAAATAGGTAACAGG
ATTATAAAAAACCAAACCTATTTTACGAGCAACCAAATTAAGCAGAAAGCAAAGAAAACAGCAATAGG
CATAAAGGAAGTTGGCCTCATGCTTGCAGCTATCCTGGCCCTCCTACTGGCTTTCTATGCCTTCTTTTAT
CTCCGACTGTCCACGAATATTGACTCTGATCTGGACCTGGATGAAGAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226002 representing NM_001171801
Red=Cloning site Green=Tags(s)

MGRKDSNTKLPVDQYRKQIGKQDYKTKPILRATKPKAEAKKTAIGIKEVGLMLAAILALLLAFYAFFY
LRLSTNIDSDLDEDED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001171801
ORF Size: 258 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001171801.1</u> , <u>NP_001165272.1</u>
RefSeq Size:	1501 bp
RefSeq ORF:	261 bp
Locus ID:	208820
UniProt ID:	<u>B2B9E1</u>
Cytogenetics:	4
MW:	9.7 kDa