

Product datasheet for **MR226003**

Triqk (NM_173746) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Triqk (NM_173746) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Triqk
Synonyms: C130086A10; Gm11818; OTTMUSG00000004461
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR226003 representing NM_173746
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGTAGAAAAGACTCTTCCAATACAAAACCTCCTGTTGATCAGTACAGAAAGCAAATAGGTAACAGG
 ATTATAAAAAACCAAACCTATTTACGAGCAACCAAATAAAAGCAGAAGCAAAGAAAACAGCAATAGG
 CATAAAGGAAGTTGGCCTCATGCTGCAGCTATCCTGGCCCTCCTACTGGCTTTCTATGCCTTCTTTAT
 CTCCGACTGTCCACGAATATTGACTCTGATCTGGACCTGGATGAAGAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226003 representing NM_173746
 Red=Cloning site Green=Tags(s)
 MGRKSSNTKLPVDQYRKQIGKQDYKTKPILRATKCLKAEAKKTAIGIKEVGLMLAAILALLLAFYAFFY
 LRLSTNIDSDLDEDED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9053_e09.zip

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_173746

ORF Size: 258 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173746.3](#), [NP_776107.2](#)


RefSeq Size: 1498 bp

RefSeq ORF: 261 bp

Locus ID: 208820

UniProt ID: [B2B9E1](#)
 Cytogenetics: 4
 MW: 9.7 kDa
 Gene Summary: May play a role in cell growth and maintenance of cell morphology.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR226003