

## Product datasheet for **MR226011L4V**

### Foxf2 (NM\_010225) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Foxf2 (NM_010225) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Foxf2
Synonyms:	Fkh20; FREAC2; LUN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010225
ORF Size:	1341 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226011).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010225.2</a> , <a href="#">NP_034355.2</a>
RefSeq Size:	2366 bp
RefSeq ORF:	1341 bp
Locus ID:	14238
UniProt ID:	<a href="#">O54743</a>
Cytogenetics:	13 A3.2
Gene Summary:	Probable transcription activator for a number of lung-specific genes (PubMed:9676429). Mediates up-regulation of the E3 ligase IRF2BPL and drives ubiquitination and degradation of CTNNB1 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)