

Product datasheet for **MR226029**

Stau1 (NM_011490) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Stau1 (NM_011490) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Stau1
Synonyms:	5830401L18Rik; AW549911; C85792; Stau
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR226029 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTATAAGCCGTGGACCTCACTCTCGGATGCAGTCCACCTACAGCTATGGCATGCGTGGAGGTGCCT
 ATCCCCCAGATACTTTTACCCATTTCCAGTCCCACCTTTACTCTACCAAGTTGAGCTCTCCGTGGGCGG
 ACAGCAGTTTAAATGGAAAGGAAAGATGAGACCACCCGTGAAACACGATGCCCTGCCCGTGCCTGAGG
 ACTCTGCAGAGTGAACCCCTGCCAGAAAGTTGGAGGTAATGGAAGAGAAGCAGAGGAAGAAAACCTCA
 ATAAATCGGAAATAAGCCAAGTGTTTGAAATTGCGCTGAAGCGGAATTTGCCTGTGAATTTGAGGTGCC
 CCGGGAGAGTGGCCACCACACATGAAGAACTTTGTGACCAGGTTTCAGTTGGGGAATTTGTAGGGGAA
 GGAGAAGGAAAAGCAAGAAGATCTCCAAGAAGATGCGGCCAGGGCTTCTGGAGCAGCTTAGGAGGC
 TGCCACCCCTCCCTGCTGTGGAGCGAGTGAAGCCAGAATCAAGAAGAAAAGTCAGCCACCTGCAAGCT
 ACAGACAGCCCCGATTATGGCCAAGGGATGAATCCTATTAGTAGACTTGCACAGATCCAGCAGGCAAAA
 AAGGAGAAGGAGCCAGAGTACATGCTCCTTACAGAACGAGGTCTCCACGTGCGCAGGAGTTTGTATGC
 AGGTAAAGGTTGGGCATCACACTGCAGAAGGAGTGGGTACCAATAAGAAGGTGGCCAAGCGTAATGCTGC
 TGAGAACATGCTGGAGATCCTGGGGTCAAAGTCCCCAGGCGCAGCCTGCCAAGCCAGCACTCAAATCA
 GAAGAGAAGACTCCAGTAAAGAAACCAGGAGACGGAAGGAAAAGTAACGTTTTTTGAACCTAGCCCTGGG
 ATGAAAATGGAAGTAGTAACAAGGACGAGGAGTTCAGGATGCCTTATCTTAGCCATCAGCAGCTGCCAGC
 TGAATTTCTCCCATGGTGCCGGAAGTTGCCAGGCTGTGGGGTTAGTCAAGGACACCACCAAAGAT
 TTCACCAGGCGAGCTCCAAATCCTGCCAAGGCAACGGTAACTGCCATGATAGCCGAGAGTTGTTGTACG
 GGGCACCTCGCCACAGCCGAGACCATTTTAAAGAGTAACATCTCTTCAGGCCAGTACCCATGGAGC
 TCGCACTAGACCCTCTGAGCAACTGTACTACCTTTCCAGAGCCAGGGATTCCAGGTTGAATACAAAGAT
 TTTCCCAAGAACAACAAGACGAGTGTGTATCTCTCATCAACTGCTCCTCACAGCCGCTCTCGTAGTC
 ATGGCATCGGCAAGGATGTGGAGTCTGTATGATATGGCTGCACTGAACATTTTAAAGCTGCTGTCTGA
 GTTGGACCAACAGAGCACAGAGATGCCAAGAACAGGAAATGGACCAGTTTCAGCGTGCGGGAGGTGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR226029 protein sequence
 Red=Cloning site Green=Tags(s)

MYKPVDPHSRMQSTYSYGMRRGAYPPRYFYFPFPVPLLYQVELSVGGQQFNGKGMPPVKHDAPARALR
 TLQSEPLPERLEVNGREAEENLNKSEISQVFEIALKRNLVNFVARESGPPHMKNFVTRVSVGEFVGE
 GEGKSKKISKKNAARAVLEQLRRLPPLPAVERVKPRIKKKQPTCKLQTAPDYQGGMNPI SRLAQIQQAK
 KEKEPEYMLLTERGLPRRREFVMQVKVGHHTAEGVGTNKKVAKRNAENMLEILGFKVPQAQPAKPKLS
 EEKTPVKKPGDGRKVTFFEPSPGDENGTSNKDEEFRMPYLSHQQLPAGILPMVPEVAQAVGVSQGHHTKD
 FTRAAPNPAKATVTAMIARELLYGGTSPTAETILKSNISSGHVPHGPRTRPSEQLYYLSRAQGFQVEYKD
 FPKNNKNECVSLINCSSQPPLVSHGIGKDVESCHDMAALNILKLLSELDQQSTEMPRTGNGPVSACGR

TRTRPLEQ**KL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

ACCN: NM_011490

ORF Size: 1470 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011490.1](#), [NM_011490.2](#), [NM_011490.3](#), [NM_011490.4](#), [NP_035620.1](#)
RefSeq Size: 3601 bp


RefSeq ORF: 1464 bp

Locus ID: 20853

UniProt ID: [Q9Z108](#)

Cytogenetics: 2 H3
MW: 54.2 kDa
Gene Summary: Binds double-stranded RNA (regardless of the sequence) and tubulin. May play a role in specific positioning of mRNAs at given sites in the cell by cross-linking cytoskeletal and RNA components, and in stimulating their translation at the site.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR226029