

Product datasheet for **MR226041L4V**

Icam5 (NM_008319) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Icam5 (NM_008319) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Icam5
Synonyms:	CD50; ICAM-3; Icam3; Tlcn; TLN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008319
ORF Size:	2751 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226041).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008319.2 , NP_032345.2
RefSeq Size:	2982 bp
RefSeq ORF:	2754 bp
Locus ID:	15898
UniProt ID:	Q60625
Cytogenetics:	9 7.7 cM
Gene Summary:	ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). [UniProtKB/Swiss-Prot Function]



[View online »](#)