

## Product datasheet for MR226057

### Eif5a (NM\_001166589) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Eif5a (NM\_001166589) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Eif5a  
**Synonyms:** AA410058; D19Wsu54; D19Wsu54e; eIF-4D; eIF-5A; eIF-5A-1; eIF-5A1; Eif4d; Eif5a1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR226057 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCAGATGATTTGGACTTCGAGACAGGAGATGCAGGGGCCACCCACCTTCCCAATGCAGTGCTCAG  
 CATTACGTAAGAATGGTTTTGTGGTGCTCAAAGGCCGCCATGTAAGATCGTCGAGATGTCTACTTCGAA  
 GACTGGCAAGCATGGCCATGCCAAGGTCCATCTGGTTGGCATTGACATTTTACTGGGAAGAAATATGAA  
 GATATCTGCCCGTCGACTCATAATATGGATGTCCCAACATCAAACGGAATGACTTCCAGCTGATTGGCA  
 TCCAGGATGGGTACCTATCCCTGCTCCAGGACAGTGGGAGGTACGAGAGGACCTTCGTCTGCCTGAAGG  
 AGACCTTGCCAAGGAGATTGAGCAGAAGTATGACTGTGGAGAAGAGATCCTGATCACAGTGTCTGCC  
 ATGACAGAGGAGGAGCTGTTGCAATCAAGCCATGGCAAAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR226057 protein sequence  
 Red=Cloning site Green=Tags(s)

MADDLDFETGDAGASATFPMQCSALRKNFVVLKGRPCKIVEMSTSKTGKHGAKVHLVGIDIFTGKKYE  
 DICPSTHNMDVPIKRNDQLIGIQDGYLSLLQDSGEVREDLRLPEGDLKEIEQKYDCGEEILITVLSA  
 MTEEAAVAIAKAMAK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** NM\_001166589

**ORF Size:** 465 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001166589.1](#), [NP\\_001160061.1](#)

**RefSeq Size:** 1455 bp

**RefSeq ORF:** 465 bp


**Locus ID:** 276770

**UniProt ID:** [P63242](#)

**Cytogenetics:** 11 B3  
**MW:** 16.8 kDa

**Gene Summary:** This gene encodes an elongation initiation factor, which participates in protein synthesis. The encoded protein also plays roles in mRNA metabolism, cell proliferation, and cell cycle control. This protein contains a modified lysine residue called hypusine, which appears to be necessary for its function. Alternatively spliced transcript variants have been described. Related pseudogenes exist on chromosomes 2, 5, and 19. [provided by RefSeq, Oct 2009]

**Product images:**



Circular map for MR226057