

Product datasheet for **MR226080**

Eif5a (NM_001166596) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Eif5a (NM_001166596) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Eif5a
Synonyms: AA410058; D19Wsu54; D19Wsu54e; eIF-4D; eIF-5A; eIF-5A-1; eIF-5A1; Eif4d; Eif5a1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR226080 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGATGATTTGGACTTCGAGACAGGAGATGCAGGGGCCACCCACCTTCCCAATGCAGTGCTCAG
CATTACGTAAGAATGGTTTTGTGGTGCTCAAAGGCCGCCATGTAAGATCGTCGAGATGTCTACTTCGAA
GACTGGCAAGCATGGCCATGCCAAGGTCCATCTGGTTGGCATTGACATTTTACTGGGAAGAAATATGAA
GATATCTGCCGTCGACTCATAATATGGATGTCCCAACATCAAACGGAATGACTCCAGCTGATTGGCA
TCCAGGATGGGTACCTATCCCTGCTCCAGGACAGTGGGGAGGTACGAGAGGACCTTCGCTGCCTGAAGG
AGACCTTGCCAAGGAGATTGAGCAGAAGTATGACTGTGGAGAAGAGATCCTGATCACAGTGTCTGCC
ATGACAGAGGAGGAGCTGTTGCAATCAAGCCATGGCAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226080 protein sequence
Red=Cloning site Green=Tags(s)

MADDLDFETGDAGASATFPMQCSALRKNFVVLKGRPCKIVEMSTSKTGKHGAKVHLVGIDIFTGKKYE
DICPSTHNMVDPNIKRNDQLIGIQDGYLSLLQDSGEVREDLRLPEGDLKEIEQKYDCGEEILITVLSA
MTEEAAVAIAKAMAK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:


ACCN: NM_001166596

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166596.1](#), [NP_001160068.1](#)

RefSeq Size: 1293 bp

RefSeq ORF: 465 bp

Locus ID: 276770


UniProt ID: [P63242](#)

Cytogenetics: 11 B3

MW: 16.8 kDa

Gene Summary: This gene encodes an elongation initiation factor, which participates in protein synthesis. The encoded protein also plays roles in mRNA metabolism, cell proliferation, and cell cycle control. This protein contains a modified lysine residue called hypusine, which appears to be necessary for its function. Alternatively spliced transcript variants have been described. Related pseudogenes exist on chromosomes 2, 5, and 19. [provided by RefSeq, Oct 2009]

Product images:



Circular map for MR226080