

Product datasheet for **MR226115L3V**

Boll (NM_001113367) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Boll (NM_001113367) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Boll
Synonyms:	4930554P13Rik; 4930597B14Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113367
ORF Size:	480 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226115).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001113367.1
RefSeq Size:	1087 bp
RefSeq ORF:	483 bp
Locus ID:	75388
UniProt ID:	Q924M5
Cytogenetics:	1 C1.2
Gene Summary:	Probable RNA-binding protein, which may be required during spermatogenesis. May act by binding to the 3' UTR of mRNAs and regulating their translation (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)