

Product datasheet for **MR226136L4V**

Ube2e3 (NM_009454) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2e3 (NM_009454) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2e3
Synonyms:	Ubce4; ubcM2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009454
ORF Size:	621 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226136).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009454.2 , NP_033480.1
RefSeq Size:	1562 bp
RefSeq ORF:	624 bp
Locus ID:	22193
UniProt ID:	P52483
Cytogenetics:	2 C3



[View online »](#)

Gene Summary:

Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-11'- and 'Lys-48'-, as well as 'Lys-63'-linked polyubiquitination (By similarity). Participates in the regulation of transepithelial sodium transport in renal cells. May be involved in cell growth arrest.[UniProtKB/Swiss-Prot Function]