

Product datasheet for **MR226146L4V**

Grik4 (NM_175481) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Grik4 (NM_175481) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Grik4
Synonyms:	6330551K01Rik; GluK4; GluRgamma1; KA-1; KA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175481
ORF Size:	2871 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226146).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_175481.5 , NP_780690.2
RefSeq Size:	3488 bp
RefSeq ORF:	2871 bp
Locus ID:	110637
UniProt ID:	Q8BMF5
Cytogenetics:	9 23.8 cM



[View online »](#)

Gene Summary:

Receptor for glutamate. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous system. The postsynaptic actions of Glu are mediated by a variety of receptors that are named according to their selective agonists. This receptor binds kainate > quisqualate > glutamate >> AMPA (By similarity).[UniProtKB/Swiss-Prot Function]