

Product datasheet for **MR226152L4V**

Dnmt3a (NM_153743) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dnmt3a (NM_153743) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dnmt3a
Synonyms:	MmullIA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153743
ORF Size:	2067 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226152).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153743.4 , NP_714965.1
RefSeq Size:	8913 bp
RefSeq ORF:	2070 bp
Locus ID:	13435
UniProt ID:	O88508
Cytogenetics:	12 A1.1



[View online »](#)

Gene Summary:

This is one of two related genes encoding de novo DNA methyltransferases, which are responsible for the establishment of DNA methylation patterns in embryos. Loss of function of this gene causes developmental defects in multiple different organ systems. There is a pseudogene for this gene located on chromosome 3. Alternatively spliced transcript variants encoding multiple isoforms have been observed. [provided by RefSeq, Nov 2012]