

Product datasheet for **MR226170L3V**

Cdc6 (NM_011799) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdc6 (NM_011799) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdc6
Synonyms:	CDC18L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011799
ORF Size:	1767 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226170).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011799.2 , NP_035929.1
RefSeq Size:	4608 bp
RefSeq ORF:	1770 bp
Locus ID:	23834
UniProt ID:	O89033
Cytogenetics:	11 D
Gene Summary:	Involved in the initiation of DNA replication. Also participates in checkpoint controls that ensure DNA replication is completed before mitosis is initiated.[UniProtKB/Swiss-Prot Function]


[View online »](#)