

Product datasheet for **MR226173L4V**

Emx2 (NM_010132) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Emx2 (NM_010132) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Emx2
Synonyms:	Pdo
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010132
ORF Size:	759 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226173).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010132.2 , NP_034262.2
RefSeq Size:	2598 bp
RefSeq ORF:	762 bp
Locus ID:	13797
UniProt ID:	Q04744
Cytogenetics:	19 56.28 cM
Gene Summary:	Transcription factor, which in cooperation with EMX2, acts to generate the boundary between the roof and archipallium in the developing brain. May function in combinations with OTX1/2 to specify cell fates in the developing central nervous system.[UniProtKB/Swiss-Prot Function]



[View online »](#)